Rigid Track

Generated by Doxygen 1.8.13

Contents

| 1 | Rigid Track Doxygen Documentation | | | | |
|---|-----------------------------------|----------|-------------|----------------------------|----|
| | 1.1 | Introdu | iction | | 1 |
| | 1.2 | Rigid 1 | rack Instal | lation | 1 |
| | 1.3 | Source | Code . | | 1 |
| 2 | File | Index | | | 3 |
| | 2.1 | File Lis | st | | 3 |
| 3 | File | Docum | entation | | 5 |
| | 3.1 | RigidTi | rack/main.c | pp File Reference | 5 |
| | | 3.1.1 | Detailed | Description | 9 |
| | | 3.1.2 | Function | Documentation | 9 |
| | | | 3.1.2.1 | calcBoardCornerPositions() | 10 |
| | | | 3.1.2.2 | calibrateGround() | 10 |
| | | | 3.1.2.3 | closeUDP() | 13 |
| | | | 3.1.2.4 | determineExposure() | 14 |
| | | | 3.1.2.5 | determineOrder() | 16 |
| | | | 3.1.2.6 | drawPositionText() | 17 |
| | | | 3.1.2.7 | getEulerAngles() | 19 |
| | | | 3.1.2.8 | loadCalibration() | 19 |
| | | | 3.1.2.9 | loadCameraPosition() | 20 |
| | | | 3.1.2.10 | loadMarkerConfig() | 21 |
| | | | 3.1.2.11 | main() | 23 |
| | | | 31212 | Mat2OPiyman() | 25 |

ii CONTENTS

| Index | | | | 73 |
|-------|---------|-------------|--------------------------|----|
| | 3.4.1 | Detailed I | Description | 71 |
| 3.4 | RigidTr | ack/RigidT | rack.h File Reference | 70 |
| | 3.3.1 | Detailed I | Description | 70 |
| 3.3 | RigidTr | ack/RigidT | rack.cpp File Reference | 70 |
| | | 3.2.3.1 | commObj | 69 |
| | 3.2.3 | Variable I | Documentation | 69 |
| | | 3.2.2.14 | testAlgorithms() | 67 |
| | | 3.2.2.13 | startTracking() | 61 |
| | | 3.2.2.12 | setReference() | 58 |
| | | 3.2.2.11 | setHeadingOffset() | 57 |
| | | 3.2.2.10 | sendDataUDP() | 56 |
| | | 3.2.2.9 | projectCoordinateFrame() | 56 |
| | | 3.2.2.8 | loadMarkerConfig() | 53 |
| | | 3.2.2.7 | loadCameraPosition() | 53 |
| | | 3.2.2.6 | loadCalibration() | 52 |
| | | 3.2.2.5 | drawPositionText() | 51 |
| | | 3.2.2.4 | determineOrder() | 49 |
| | | 3.2.2.3 | determineExposure() | 47 |
| | | 3.2.2.2 | closeUDP() | 46 |
| | | 3.2.2.1 | calibrateGround() | 44 |
| | 3.2.2 | Function | Documentation | 43 |
| | 3.2.1 | Detailed I | Description | 43 |
| 3.2 | RigidTr | rack/main.h | File Reference | 41 |
| | | 3.1.3.2 | Rmat | 40 |
| | | 3.1.3.1 | commObj | 40 |
| | 3.1.3 | Variable I | Documentation | 40 |
| | | 3.1.2.18 | testAlgorithms() | 37 |
| | | 3.1.2.17 | startTracking() | 31 |
| | | 3.1.2.16 | setReference() | 28 |
| | | 3.1.2.15 | setHeadingOffset() | 28 |
| | | 3.1.2.14 | sendDataUDP() | 27 |
| | | 3.1.2.13 | projectCoordinateFrame() | 26 |

Chapter 1

Rigid Track Doxygen Documentation

1.1 Introduction

Rigid Track is a software that provides, combined with an OptiTrack camera, the pose estimation of one object in three dimensional space. This is achieved with only one camera in combination with reflective markers. Those are attached to the object ought to be tracked. The accuracy in the range of millimeters and the high update rate of 100 Hz enable use cases for fast and agile objects. The main application is navigation for drones that rely on high precision position data. Where GPS is not available, e.g. indoors or due to a lacking GPS receiver, this setup substitutes for it. Another use case is the pure pose logging when the drone does not depend on the position, e.g. when it is remote piloted by hand. While this setup contains one OptiTrack Flex 3 camera, every other model of OptiTrack should work, despite not tested. With better camera models, e.g. the Prime Series, even outdoor usage is possible. When the capabilities are not sufficient please refer to OptiTracks Software Motive. But keep in mind that this solution needs at least 3 cameras as Rigid Track works with only one.

1.2 Rigid Track Installation

Start the RigidTrack_setup.exe from the enclosed SD card and follow the instructions given in the installation assistant. Default parameters like installation directory or shortcuts to be created can be chosen. But normally clicking Next and keeping the default values should be sufficient. When the installation is completed a shortcut in the start menu and the desktop can be used to start Rigid Track. The program is then successfully installed in C:/Program Files (x86)/TU Munich FSD/Rigid Track.

1.3 Source Code

The most interesting file for you is main.cpp. It contains the relevant functions for pose estimation. Camera calibration and other functional aspects are also implemented there. The GUI program code is found in RigidTrack.cpp. communication.cpp deals only with communication from main.cpp to the GUI.

Chapter 2

File Index

2.1 File List

Here is a list of all documented files with brief descriptions:

| RigidTrack/_modelest.h |
|---|
| RigidTrack/communication.cpp |
| RigidTrack/communication.h |
| RigidTrack/main.cpp |
| Rigid Track main file that contains most functionallity |
| RigidTrack/main.h |
| Header file for main.cpp |
| RigidTrack/precomp.hpp ?? |
| RigidTrack/resource.h?? |
| RigidTrack/RigidTrack.cpp |
| Rigid Track GUI source that contains functions for GUI events |
| RigidTrack/RigidTrack.h |
| Rigid Track GUI source header with Qt Signals and Slots |
| RigidTrack/supportcode.cpp |
| RigidTrack/supportcode.h |

File Index

Chapter 3

File Documentation

3.1 RigidTrack/main.cpp File Reference

Rigid Track main file that contains most functionallity.

```
#include "RigidTrack.h"
#include "main.h"
#include "communication.h"
#include "cameralibrary.h"
#include "modulevector.h"
#include "modulevectorprocessing.h"
#include "coremath.h"
#include <QtWidgets/QApplication>
#include <QDesktopServices>
#include <QInputDialog>
#include <QUrl>
#include <QThread>
#include <QUdpSocket>
#include <QFileDialog>
#include <opencv\cv.h>
#include "opencv2\core.hpp"
#include "opencv2\calib3d.hpp"
#include <opencv2/imgproc/imgproc.hpp>
#include <opencv2/calib3d/calib3d.hpp>
#include <opencv2/highqui/highqui.hpp>
#include <opencv2\video\tracking.hpp>
#include <fstream>
#include <windows.h>
#include <conio.h>
#include <tchar.h>
#include <stdio.h>
#include <iostream>
#include <stdarg.h>
#include <ctype.h>
#include <stdlib.h>
#include <gl/glu.h>
#include <sstream>
#include <time.h>
#include <cmath>
#include <vector>
```

```
#include <algorithm>
#include <random>
#include <thread>
#include <strsafe.h>
```

Include dependency graph for main.cpp:



Functions

• int main (int argc, char *argv[])

main initialises the GUI and values for the marker position etc

- QPixmap Mat2QPixmap (cv::Mat src)
- void calcBoardCornerPositions (Size boardSize, float squareSize, std::vector< Point3f > &corners)
- void getEulerAngles (Mat &rotCamerMatrix, Vec3d &eulerAngles)
- int startTracking ()
- · void startStopCamera ()

Start or stop the tracking depending on if the camera is currently running or not.

- int setReference ()
- int calibrateCamera ()

Start the camera calibration routine that computes the camera matrix and distortion coefficients.

- void loadCalibration (int method)
- void testAlgorithms ()
- void projectCoordinateFrame (Mat pictureFrame)
- · void setUpUDP ()

Open the UDP ports for communication.

- void setHeadingOffset (double d)
- void sendDataUDP (cv::Vec3d &Position, cv::Vec3d &Euler)
- void closeUDP ()
- void loadMarkerConfig (int method)
- void drawPositionText (cv::Mat &Picture, cv::Vec3d &Position, cv::Vec3d &Euler, double error)
- void loadCameraPosition ()
- int determineExposure ()
- void determineOrder ()
- int calibrateGround ()

Variables

commObject commObj

class that handles the communication from main.cpp to the GUI

• bool safetyEnable = false

is the safety feature enabled

• bool safety2Enable = false

is the second receiver enabled

double safetyBoxLength = 1.5

length of the safety area cube in meters

• int safetyAngle = 30

bank and pitch angle protection in degrees

• bool exitRequested = true

variable if tracking loop should be exited

int invertZ = 1

dummy variable to invert Z direction on request

• double frameTime = 0.01

100 Hz CoSy rate, is later on replaced with the hardware timestamp delivered by the camera

• double timeOld = 0.0

old time for finite differences velocity calculation. Is later on replaced with the hardware timestamp delivered by the camera

double timeFirstFrame = 0

Time stamp of the first frame. This value is then subtracted for every other frame so the time in the log start at zero.

Vec3d position = Vec3d()

position vector x,y,z for object position in O-CoSy, unit is meter

Vec3d eulerAngles = Vec3d()

Roll Pitch Heading in this order, units in degrees.

Vec3d positionOld = Vec3d()

old position in O-CoSy for finite differences velocity calculation

Vec3d velocity = Vec3d()

velocity vector of object in o-CoSy in respect to o-CoSy

Vec3d posRef = Vec3d()

initial position of object in camera CoSy

Vec3d eulerRef = Vec3d()

initial euler angle of object respectivley to camera CoSy

• double headingOffset = 0

heading offset variable for aligning INS heading with tracking heading

• int intlntensity = 15

max infrared spot light intensity is 15 1-6 is strobe 7-15 is continuous 13 and 14 are meaningless

• int intExposure = 1

max is 480 increase if markers are badly visible but should be determined automatically during setReference()

• int intFrameRate = 100

CoSy rate of camera, maximum is 100 fps.

• int intThreshold = 200

threshold value for marker detection. If markers are badly visible lower this value but should not be necessary

Mat Rmat = (cv::Mat <double>(3, 1) << 0.0, 0.0, 0.0)

Rotation, translation etc. matrix for PnP results.

Mat RmatRef = (cv::Mat_<double>(3, 3) << 1., 0., 0., 0., 1., 0., 0., 1.)

reference rotation matrix from camera CoSy to marker CoSy

Mat M CN = cv::Mat <double>(3, 3)

rotation matrix from camera to ground, fixed for given camera position

Mat M_HeadingOffset = cv::Mat_<double>(3, 3)

rotation matrix that turns the ground system to the INS magnetic heading for alignment

Mat Rvec = (cv::Mat_<double>(3, 1) << 0.0, 0.0, 0.0)

rotation vector (axis-angle notation) from camera CoSy to marker CoSy

• Mat Tvec = (cv::Mat_<double>(3, 1) << 0.0, 0.0, 0.0)

translation vector from camera CoSy to marker CoSy in camera CoSy

Mat RvecOriginal

initial values as start values for algorithms and algorithm tests

Mat TvecOriginal

initial values as start values for algorithms and algorithm tests

• bool useGuess = true

set to true and the algorithm uses the last result as starting value

• int methodPNP = 0

solvePNP algorithm 0 = iterative 1 = EPNP 2 = P3P 4 = UPNP //!< 4 and 1 are the same and not implemented correctly by OpenCV

• int numberMarkers = 4

number of markers. Is loaded during start up from the marker configuration file

std::vector< Point3d > list_points3d

marker positions in marker CoSy

std::vector< Point2d > list_points2d

marker positions projected in 2D in camera image CoSy

std::vector< Point2d > list_points2dOld

marker positions in previous picture in 2D in camera image CoSy

std::vector< double > list_points2dDifference

difference of the old and new 2D marker position to determine the order of the points

std::vector< Point2d > list_points2dProjected

3D marker points projected to 2D in camera image CoSy with the algorithm projectPoints

std::vector< Point2d > list_points2dUnsorted

marker points in 2D camera image CoSy, sorted with increasing x (camera image CoSy) but not sorted to correspond with list_points3d

std::vector < Point3d > coordinateFrame

coordinate visualisazion of marker CoSy

std::vector< Point2d > coordinateFrameProjected

marker CoSy projected from 3D to 2D camera image CoSy

int pointOrderIndices [] = { 0, 1, 2, 3 }

old correspondence from list_points3d and list_points_2d

int pointOrderIndicesNew [] = { 0, 1, 2, 3 }

new correspondence from list_points3d and list_points_2d

double currentPointDistance = 5000

distance from the projected 3D points (hence in 2d) to the real 2d marker positions in camera image CoSy

• double minPointDistance = 5000

minimum distance from the projected 3D points (hence in 2d) to the real 2d marker positions in camera image CoSy

• int currentMinIndex = 0

helper variable set to the point order that holds the current minimum point distance

bool gotOrder = false

order of the list_points3d and list_points3d already tetermined or not, has to be done once

bool camera_started = false

variable thats needed to exit the main while loop

· Mat cameraMatrix

camera matrix of the camera

Mat distCoeffs

distortion coefficients of the camera

· Core::DistortionModel distModel

distortion model of the camera

QUdpSocket * udpSocketObject

socket for the communication with receiver 1

QUdpSocket * udpSocketSafety

socket for the communication with safety receiver

QUdpSocket * udpSocketSafety2

socket for the communication with receiver 3

QHostAddress IPAdressObject = QHostAddress("127.0.0.1")

IPv4 adress of receiver 1.

QHostAddress IPAdressSafety = QHostAddress("192.168.4.1")

IPv4 adress of safety receiver.

QHostAddress IPAdressSafety2 = QHostAddress("192.168.4.4")

IPv4 adress of receiver 2.

• int portObject = 9155

Port of receiver 1.

int portSafety = 9155

Port of the safety receiver.

• int portSafety2 = 9155

Port of receiver 2.

QByteArray datagram

data package that is sent to receiver 1 and 2

QByteArray data

data package that's sent to the safety receiver

• const int BACKBUFFER_BITSPERPIXEL = 8

8 bit per pixel and greyscale image from camera

· std::string strBuf

buffer that holds the strings that are sent to the Qt GUI

• std::stringstream ss

stream that sends the strBuf buffer to the Qt GUI

QString logFileName

Filename for the logfiles.

std::string logName

Filename for the logfiles as standard string.

SYSTEMTIME logDate

Systemtime struct that saves the current date and time thats needed for the log file name creation.

std::ofstream logfile

file handler for writing the log file

3.1.1 Detailed Description

Rigid Track main file that contains most functionallity.

This file contains allmost all functional code for pose estimation, calibration and so on. The GUI related part is in RigidTrack.cpp and the communication from main.cpp to GUI is done with the commObj class from communication.cpp.

Author

Florian J.T. Wachter

Version

1.0

Date

April, 8th 2017

3.1.2 Function Documentation

3.1.2.1 calcBoardCornerPositions()


Calculate the chess board corner positions, used for the camera calibration.

Parameters

| in | boardSize | denotes how many squares are in each direction. |
|-----|------------|---|
| in | squareSize | is the square length in millimeters. |
| out | corners | returns the square corners in millimeters. |

Definition at line 229 of file main.cpp.

Here is the caller graph for this function:



3.1.2.2 calibrateGround()

```
int calibrateGround ( )
```

Get the pose of the camera w.r.t the ground calibration frame. This frame sets the navigation frame for later results. The pose is averaged over 200 samples and then saved in the file referenceData.xml. This routine is basically the same as setReference.


Definition at line 1563 of file main.cpp.

```
1564 {
1565
          //! initialize the variables with starting values
1566
         gotOrder = false;
         posRef = 0;
1567
1568
         eulerRef = 0:
         RmatRef = 0;
1569
1570
         Rvec = RvecOriginal;
1571
         Tvec = TvecOriginal;
1572
1573
         determineExposure();
1574
1575
         ss.str("");
1576
         commObj.addLog("Started ground calibration");
1577
1578
         CameraLibrary_EnableDevelopment();
1579
          //! Initialize Camera SDK
1580
         CameraLibrary::CameraManager::X();
1581
1582
         //! At this point the Camera SDK is actively looking for all connected cameras and will initialize
1583
         //! them on it's own.
1584
1585
          //! Get a connected camera =========
1586
         CameraManager::X().WaitForInitialization();
1587
         Camera *camera = CameraManager::X().GetCamera();
1588
1589
         //! If no device connected, pop a message box and exit ==--
1590
         if (camera == 0)
1591
1592
              commObj.addLog("No camera found!");
1593
              return 1;
1594
1595
1596
         //! Determine camera resolution to size application window ==----
         int cameraWidth = camera->Width();
int cameraHeight = camera->Height();
1597
1598
         camera->GetDistortionModel(distModel);
1599
1600
         cv::Mat matFrame(cv::Size(cameraWidth, cameraHeight), CV_8UC1);
1601
1602
          //! Set camera mode to precision mode, it directly provides marker coordinates
1603
         camera->SetVideoType(Core::PrecisionMode);
1604
1605
          //! Start camera output ==--
1606
         camera->Start();
1607
1608
           /! Turn on some overlay text so it's clear things are
1609
         //! working even if there is nothing in the camera's view. ===---
1610
          //! Set some other parameters as well of the camera
1611
         camera->SetTextOverlay(true);
1612
         camera->SetFrameRate(intFrameRate);
1613
         camera->SetIntensity(intIntensity);
1614
         camera->SetIRFilter(true);
1615
         camera->SetContinuousIR(false);
1616
         camera->SetHighPowerMode(false);
1617
         //! sample some frames and calculate the position and attitude. then average those values and use that
1618
       as zero position
1619
        int numberSamples = 0;
1620
          int numberToSample = 200;
1621
         double projectionError = 0;
1622
1623
         while (numberSamples < numberToSample)</pre>
1624
         {
1625
              //! Fetch a new frame from the camera ===---
             Frame *frame = camera->GetFrame();
1626
1627
1628
              if (frame)
1629
                  //! Ok, we've received a new frame, lets do something
1630
1631
                  //! with it.
1632
                  if (frame->ObjectCount() == numberMarkers)
1633
1634
                      //!for(int i=0; i<frame->ObjectCount(); i++)
1635
                      for (int i = 0; i < numberMarkers; i++)</pre>
1636
                          cObject *obj = frame->Object(i);
1637
                          list_points2dUnsorted[i] = cv::Point2d(obj->X(), obj->Y());
1638
1639
1640
1641
                      if (gotOrder == false)
1642
1643
                          determineOrder();
1644
1645
1646
                      //! sort the 2d points with the correct indices as found in the preceeding order
       {\tt determination}\ {\tt algorithm}
                      for (int w = 0; w < numberMarkers; w++)</pre>
1647
1648
```

```
1649
                          list_points2d[w] = list_points2dUnsorted[
      pointOrderIndices[w]];
1650
1651
                      list_points2dOld = list_points2dUnsorted;
1652
1653
                      //!Compute the pose from the 3D-2D corresponses
                      solvePnP(list_points3d, list_points2d,
1654
      cameraMatrix, distCoeffs, Rvec, Tvec, useGuess,
      methodPNP);
1655
                      //! project the marker 3d points with the solution into the camera image CoSy and calculate
1656
       difference to true camera image
1657
                     projectPoints(list_points3d, Rvec, Tvec,
      cameraMatrix, distCoeffs, list_points2dProjected);
1658
                     projectionError = norm(list_points2dProjected,
      list_points2d);
1659
1660
                      if (projectionError > 3)
1661
1662
                          commObj.addLog("Reprojection error is bigger than 3 pixel. Correct marker
       configuration loaded?\nMarker position measured precisely?");
1663
                         frame->Release();
1664
                          return 1;
1665
1666
1667
                      double maxValue = 0;
1668
                      double minValue = 0;
1669
                      minMaxLoc(Tvec.at<double>(2), &minValue, &maxValue);
1670
1671
                      if (maxValue > 10000 || minValue < 0)</pre>
1672
1673
1674
1675
                          commObj.addLog("Negative z distance, thats not possible. Start the set zero
       routine again and check marker configurations.");
1676
                         frame->Release();
1677
                          return 1;
1678
1679
                     if (norm(positionOld) - norm(Tvec) < 0.05) //!<Iterative Method needs time</pre>
1680
       to converge to solution
1681
                     {
                          add(posRef, Tvec, posRef);
1682
                          add(eulerRef, Rvec, eulerRef); //! That are not the values of yaw,
1683
       roll and pitch yet! Rodriguez has to be called first.
1684
                          numberSamples++;
                                             //!<-- one sample more :D
1685
                          commObj.progressUpdate(numberSamples * 100 / numberToSample);
1686
                     positionOld = Tvec;
1687
1688
1689
                      Mat cFrame(480, 640, CV_8UC3, Scalar(0, 0, 0));
1690
                      for (int i = 0; i < numberMarkers; i++)</pre>
1691
1692
                          circle(cFrame, Point(list_points2d[i].x,
      list_points2d[i].y), 6, Scalar(0, 225, 0), 3);
1693
1694
                      projectCoordinateFrame(cFrame);
1695
                      projectPoints(list_points3d, Rvec, Tvec,
      cameraMatrix, distCoeffs, list_points2d);
1696
                      for (int i = 0; i < numberMarkers; i++)</pre>
1697
                          circle(cFrame, Point(list_points2d[i].x,
1698
      list_points2d[i].y), 3, Scalar(225, 0, 0), 3);
1699
1700
1701
                      QPixmap QPFrame;
1702
                      QPFrame = Mat2QPixmap(cFrame);
1703
                      commObj.changeImage(OPFrame);
1704
                      OCoreApplication::processEvents();
1705
1706
1707
                  frame->Release();
1708
             }
1709
         //! Release camera ==--
1710
1711
         camera->Release();
1712
1713
         \ensuremath{//!}\ensuremath{\text{Divide}} by the number of samples to get the mean of the reference position
1714
         divide(posRef, numberToSample, posRef);
         divide(eulerRef, numberToSample, eulerRef); //!< eulerRef is here in Axis Angle
1715
       notation
1716
1717
         Rodrigues(eulerRef, RmatRef);
                                                        //!< axis angle to rotation matrix
1718
1719
         getEulerAngles(RmatRef, eulerRef); //!< rotation matrix to euler</pre>
         ss.str("");
1720
1721
         ss << "RmatRef is:\n";
```

```
ss << RmatRef << "\n";
             ss << "Reference Position is:\n";
ss << posRef << "[mm] \n";
ss << "Reference Euler angles are:\n";
ss << eulerRef << "[deg] \n";
1723
1724
1725
1726
1727
         //! Save the obtained calibration coefficients in a file for later use
  QString fileName = QFileDialog::getSaveFileName(nullptr, "Save ground calibration file", "
referenceData.xml", "Calibration File (*.xml);; All Files (*)");
1728
1729
              FileStorage fs(fileName.toUtf8().constData(), FileStorage::WRITE);
fs << "M_NC" << RmatRef;
fs << "eulerRef" << eulerRef;
1730
1731
1732
              strBuf = fs.releaseAndGetString();
1733
1734
              commObj.changeStatus(QString::fromStdString(strBuf));
1735
              commObj.addLog("Saved ground calibration!");
1736
1737
              commObj.progressUpdate(0);
              return 0;
1738 }
```

Here is the call graph for this function:



3.1.2.3 closeUDP()


```
void closeUDP ( )
```

Close the UDP ports again to release network interfaces etc. If this is not done the network resources are still occupied and the program can't exit properly.

Definition at line 1173 of file main.cpp.

```
1181
         if (udpSocketSafety->isOpen())
1182
1183
             udpSocketSafety->close();
1184
        }
1185
         if (udpSocketSafety2->isOpen())
1186
1187
        {
1188
             udpSocketSafety2->close();
1189
         commObj.addLog("Closed all UDP ports.");
1190
1191 }
```

Here is the caller graph for this function:



3.1.2.4 determineExposure()

```
int determineExposure ( )
```

Get the optimal exposure for the camera. For that find the minimum and maximum exposure were the right number of markers are detected. Then the mean of those two values is used as exposure.


Definition at line 1362 of file main.cpp.

```
1363 {
1364
         //! For OptiTrack Ethernet cameras, it's important to enable development mode if you
1365
         //! want to stop execution for an extended time while debugging without disconnecting
1366
        //! the Ethernet devices. Lets do that now:
1367
1368
        CameraLibrary_EnableDevelopment();
1369
1370
         //! Initialize Camera SDK ==-
1371
        CameraLibrary::CameraManager::X();
1372
1373
         //! At this point the Camera SDK is actively looking for all connected cameras and will initialize
1374
        //! them on it's own.
1375
1376
         //! Get a connected camera =======
1377
        CameraManager::X().WaitForInitialization();
1378
        Camera *camera = CameraManager::X().GetCamera();
1379
1380
         //! If no device connected, pop a message box and exit ==--
1381
        if (camera == 0)
1382
        {
1383
             commObj.addLog("No camera found!");
1384
1385
1386
        //! Determine camera resolution to size application window ==----
1387
1388
        int cameraWidth = camera->Width();
        int cameraHeight = camera->Height();
1389
1390
1391
        camera->SetVideoType(Core::PrecisionMode); //! set the camera mode to precision mode, it used
      greyscale imformation for marker property calculations
1392
1393
                                                     //! Start camera output ==--
1394
        camera->Start();
1395
```

```
1396
          //! Turn on some overlay text so it's clear things are
          //! working even if there is nothing in the camera's view. ===---
1397
1398
          camera->SetTextOverlay(true);
1399
          camera->SetExposure(intExposure);
                                                    //! set the camera exposure
          camera->SetIntensity(intIntensity); //! set the camera infrared LED intensity camera->SetFrameRate(intFrameRate); //! set the camera framerate to 100 Hz
1400
1401
1402
          camera->SetIRFilter(true); //! enable the filter that blocks visible light and only passes infrared
         camera->SetHighPowerMode(true); //! enable high power mode of the leds
camera->SetContinuousIR(false); //! enable continuous LED light
1403
1404
          camera->SetThreshold(intThreshold); //! set threshold for marker detection
1405
1406
1407
          //!set exposure such that num markers are visible
          int numberObjects = 0; //! Number of objects (markers) found in the current picture with the given
1408
        exposure
          int minExposure = 1; //! exposure when objects detected the first time is numberMarkers int maxExposure = 480; //! exposure when objects detected is first time numberMarkers+1
1409
1410
          intExposure = minExposure; //! set the exposure to the smallest value possible
1411
                                    //! if the markers arent found after numberTries then there might be no markers
1412
          int numberTries = 0;
       at all in the real world
1413
1414
                                     //! Determine minimum exposure, hence when are numberMarkers objects detected
          camera->SetExposure(intExposure);
1415
1416
          while (numberObjects != numberMarkers && numberTries < 48)</pre>
1417
1418
               //! get a new camera frame
              Frame *frame = camera->GetFrame();
1419
              if (frame) //! frame received
1420
1421
                  numberObjects = frame->ObjectCount();    //! how many objects are detected in the image
if (numberObjects == numberMarkers) { minExposure =
1422
1423
      intExposure; frame->Release(); break; } //! if the right amount if markers is found, exit while
1424
                   //! not the right amount of markers was found so increase the exposure and try again
                   numberTries++;
intExposure += 10;
1425
1426
                   camera->SetExposure(intExposure);
1427
                   ss.str("");
1428
                   ss << "Exposure: " << intExposure << "\t";
ss << "Objects found: " << numberObjects;
1429
1430
1431
                   commObj.addLog(QString::fromStdString(ss.str()));
1432
                   frame->Release();
              }
1433
1434
         }
1435
1436
          //! Now determine maximum exposure, hence when are numberMarkers+1 objects detected
1437
         numberTries = 0;
                               //! if the markers arent found after numberTries then there might be no markers at
       all in the real world \,
1438
          intExposure = maxExposure;
          camera->SetExposure(intExposure);
1439
1440
          numberObjects = 0;
1441
          while (numberObjects != numberMarkers && numberTries < 48)</pre>
1442
1443
              Frame *frame = camera->GetFrame();
1444
              if (frame)
1445
1446
                   numberObjects = frame->ObjectCount(); //! how many objects are detected in the image
                   if (numberObjects == numberMarkers) { maxExposure
1447
      intExposure; frame->Release(); break; } //! if the right amount if markers is found, exit while
1448
1449
                   //! not the right amount of markers was found so decrease the exposure and try again
1450
                   intExposure -= 10;
1451
                   numberTries++;
1452
                   camera->SetExposure(intExposure);
                   ss.str("");
1453
                   ss << "Exposure: " << intExposure << "\t";
ss << "Objects found: " << numberObjects;
1454
1455
1456
                   commObj.addLog(QString::fromStdString(ss.str()));
1457
                   frame->Release();
1458
1459
1460
          //! set the exposure to the mean of min and max exposure determined
1461
1462
          camera->SetExposure((minExposure + maxExposure) / 2.0);
1463
1464
          //! and now check if the correct amount of markers is detected with that new value
1465
          while (1)
1466
1467
              Frame *frame = camera->GetFrame():
1468
               if (frame)
1469
1470
                   numberObjects = frame->ObjectCount(); //! how many objects are detected in the image
1471
                   if (numberObjects != numberMarkers) //! are all markers and not more or less
       detected in the image
1472
                  {
1473
                       frame->Release();
```

```
commObj.addLog("Was not able to detect the right amount of markers.");
1475
                       //! Release camera =
1476
                       camera->Release();
1477
                       return 1;
1478
1479
                   else //! all markers and not more or less are found
1480
1481
                       frame->Release();
1482
                       intExposure = (minExposure + maxExposure) / 2.0;
                       commObj.addLog("Found the correct number of markers.");
commObj.addLog("Exposure set to:");
1483
1484
1485
                       commObj.addLog(QString::number(intExposure));
1486
                       break;
1487
1488
1489
1490
          camera->Release();
1491
1492
          return 0;
1493
1494 }
```

Here is the caller graph for this function:



3.1.2.5 determineOrder()

```
void determineOrder ( )
```


Compute the order of the marker points in 2D so they are the same as in the 3D array. Hence marker 1 must be in first place for both, list_points2d and list_points3d.

Definition at line 1498 of file main.cpp.

```
1499 {
1500
         //! determine the 3D-2D correspondences that are crucial for the PnP algorithm
1501
         //! \  \, {\rm Try \ every \ possible \ correspondence \ and \ solve \ PnP}
         //! Then project the 3D marker points into the 2D camera image and check the difference
1502
1503
         //! between projected points and points as seen by the camera
1504
         //! the corresponce with the smallest difference is probably the correct one
1505
1506
              //! the difference between true 2D points and projected points is super big
         minPointDistance = 5000;
std::sort(pointOrderIndices, pointOrderIndices + 4);
1507
1508
1509
1510
         //! now try every possible permutation of correspondence
1511
              //! reset the starting values for solvePnP
1512
1513
              Rvec = RvecOriginal;
1514
              Tvec = TvecOriginal;
1515
1516
              //! sort the 2d points with the current permutation
```

```
1517
              for (int m = 0; m < numberMarkers; m++)</pre>
1518
1519
                 list_points2d[m] = list_points2dUnsorted[
      pointOrderIndices[m]];
1520
             }
1521
             //! Call solve PNP with P3P since its more robust and sufficient for start value determination
1522
1523
             solvePnP(list_points3d, list_points2d,
      cameraMatrix, distCoeffs, Rvec, Tvec, useGuess, SOLVEPNP_P3P);
1524
             //! set the current difference of all point correspondences to zero
1525
1526
             currentPointDistance = 0;
1527
1528
             //! project the 3D points with the solvePnP solution onto 2D
1529
             projectPoints(list_points3d, Rvec, Tvec,
      cameraMatrix, distCoeffs, list_points2dProjected);
1530
1531
              //! now compute the absolute difference (error)
             for (int n = 0; n < numberMarkers; n++)</pre>
1532
                  currentPointDistance += norm(list_points2d[n] -
      list_points2dProjected[n]);
1535
             }
1536
1537
             //! if the difference with the current permutation is smaller than the smallest value till now
             //! it is probably the more correct permutation
1538
1539
              if (currentPointDistance < minPointDistance)</pre>
1540
1541
                  minPointDistance = currentPointDistance;
                                                              //!< set the
       smallest value of difference to the current one
                 for (int b = 0; b < numberMarkers; b++)</pre>
1542
                                                             //!< now safe the better permutation
1543
                  {
1544
                      pointOrderIndicesNew[b] = pointOrderIndices[b];
1545
1546
             }
1547
1548
1549
1550
         //! try every permutation
1551
         while (std::next_permutation(pointOrderIndices,
      pointOrderIndices + 4));
1552
         //! now that the correct order is found assign it to the indices array for (int w = 0; w < numberMarkers; w++)
1553
1554
1555
1556
             pointOrderIndices[w] = pointOrderIndicesNew[w];
1557
1558
         gotOrder = true;
1559 }
```

Here is the caller graph for this function:



3.1.2.6 drawPositionText()

```
cv::Vec3d & Position,
cv::Vec3d & Euler,
double error )
```

Draw the position, attitude and reprojection error in the picture.


Parameters

| in | Picture | is the camera image in OpenCV matrix format. |
|----|---|---|
| in | Position | is the position of the tracked object in navigation CoSy. |
| in | Euler are the Euler angles with respect to the navigation frame | |
| in | error | is the reprojection error of the pose estimation. |

Definition at line 1315 of file main.cpp.

```
1316 {
1317
          ss.str("");
ss << "X: " << Position[0] << " m";</pre>
1318
          putText(Picture, ss.str(), cv::Point(200, 440), 1, 1, cv::Scalar(255, 255, 255));
1320
          ss.str("");
ss << "Y: " << Position[1] << " m";</pre>
1321
1322
          putText(Picture, ss.str(), cv::Point(200, 455), 1, 1, cv::Scalar(255, 255, 255));
1323
1324
1325
          ss.str("");
ss << "Z: " << Position[2] << " m";</pre>
1326
1327
          putText(Picture, ss.str(), cv::Point(200, 470), 1, 1, cv::Scalar(255, 255, 255));
1328
1329
          ss.str("");
          ss << "Heading: " << Euler[2] << " deg";
1330
1331
          putText(Picture, ss.str(), cv::Point(350, 440), 1, 1, cv::Scalar(255, 255, 255));
1332
          ss.str("");
ss << "Pitch: " << Euler[1] << " deg";</pre>
1333
1334
          putText(Picture, ss.str(), cv::Point(350, 455), 1, 1, cv::Scalar(255, 255, 255));
1335
1336
1337
1338
          ss << "Roll: " << Euler[0] << " deg";
          putText(Picture, ss.str(), cv::Point(350, 470), 1, 1, cv::Scalar(255, 255, 255));
1339
1340
          ss.str("");
ss << "Error: " << error << " px";</pre>
1341
1342
          putText(Picture, ss.str(), cv::Point(10, 470), 1, 1, cv::Scalar(255, 255, 255));
1343
1344 }
```

Here is the caller graph for this function:



3.1.2.7 getEulerAngles()

Get the euler angles from a rotation matrix


Parameters

| in | rotCamerMatrix | is a projection matrix, here normally only the extrinsic values. |
|-----|----------------|--|
| out | eulerAngles | contains the Euler angles that result in the same rotation matrix as rotCamerMatrix. |

Definition at line 241 of file main.cpp.

```
241
242
243
         Mat cameraMatrix, rotMatrix, transVect, rotMatrixX, rotMatrixY, rotMatrixZ;
244
         double* _r = rotCamerMatrix.ptr<double>();
         double projMatrix[12] = { _r[0],_r[1],_r[2],0,
    _r[3],_r[4],_r[5],0,
    _r[6],_r[7],_r[8],0 };
245
246
247
248
249
         decomposeProjectionMatrix(Mat(3, 4, CV_64FC1, projMatrix),
250
           cameraMatrix,
251
              rotMatrix,
252
             transVect,
253
             rotMatrixX,
rotMatrixY,
254
255
             rotMatrixZ,
             eulerAngles);
257 }
```

Here is the caller graph for this function:



3.1.2.8 loadCalibration()

Load a previously saved camera calibration from a file.

Parameters

in method whether or not load the camera calibration from calibration.xml. If ==0 then yes, if != 0 then let the user select a different file.

Definition at line 923 of file main.cpp.

```
923
                                              {
924
925
         QString fileName;
926
         if (method == 0)
927
928
               fileName = "calibration.xml";
929
930
         else
931
         {
932
              fileName = QFileDialog::getOpenFileName(nullptr, "Choose a previous saved calibration file", "", "
       Calibration Files (*.xml);; All Files (*)");
933
               if (fileName.length() == 0)
934
                   fileName = "calibration.xml";
935
936
              }
937
938
         FileStorage fs;
939
         fs.open(fileName.toUtf8().constData(), FileStorage::READ);
         fs["CameraMatrix"] >> cameraMatrix;
fs["DistCoeff"] >> distCoeffs;
commobj.addLog("Loaded calibration from file:");
940
941
942
         commObj.addLog(fileName);
ss.str("");
943
944
         ss << "\nCamera Matrix is" << "\n" << cameraMatrix << "\n";
ss << "\nDistortion Coefficients are" << "\n" << distCoeffs << "\n";
945
946
947
         commObj.addLog(QString::fromStdString(ss.str()));
948 }
```

Here is the caller graph for this function:



3.1.2.9 loadCameraPosition()

```
void loadCameraPosition ( )
```

Load the rotation matrix from camera CoSy to ground CoSy It is determined during calibrateGround() and stays the same once the camera is mounted and fixed.

Definition at line 1348 of file main.cpp.

```
1349 {
1350
           //! Open the referenceData.xml that contains the rotation from camera CoSy to ground CoSy
1351
           FileStorage fs;
1352
           fs.open("referenceData.xml", FileStorage::READ);
          fs.Upen( referenceData.
fs["M_NC"] >> M_CN;
fs["M_NC"] >> RmatRef;
fs["posRef"] >> posRef;
1353
1354
1355
1356
           fs["eulerRef"] >> eulerRef;
1357
           commObj.addLog("Loaded reference pose.");
1358 }
```

Here is the caller graph for this function:



3.1.2.10 loadMarkerConfig()

Load a marker configuration from file. This file has to be created by hand, use the standard marker configuration file as template.

Parameters

| in | method | whether or not load the configuration from the markerStandard.xml. If ==0 load it, if != 0 let the | |
|----|--------|--|--|
| | | user select a different file. | |


Definition at line 1195 of file main.cpp.

```
1196 {
         QString fileName;
1197
1198
         //! during start up of the programm load the standard marker configuration
1199
         if (method == 0)
1200
1201
              //! open the standard marker configuration file
1202
             FileStorage fs;
             fs.open("markerStandard.xml", FileStorage::READ);
1203
1204
1205
              //! copy the values to the respective variables
1206
             fs["numberMarkers"] >> numberMarkers;
1207
1208
              //! inizialise vectors with correct length depending on the number of markers
             list_points3d = std::vector<Point3d> (numberMarkers);
list_points2d = std::vector<Point2d> (numberMarkers);
1209
1210
              list_points2dOld = std::vector<Point2d>(numberMarkers);
1211
1212
              list_points2dDifference = std::vector<double>(
      numberMarkers);
1213
             list_points2dProjected = std::vector<Point2d>(
      numberMarkers);
1214
              list_points2dUnsorted = std::vector<Point2d>(
      numberMarkers);
1215
```

```
//! save the marker locations in the points3d vector
              fs["list_points3d"] >> list_points3d;
1217
1218
              fs.release();
              commObj.addLog("Loaded marker configuration from file:");
1219
              commObj.addLog(fileName);
1220
1221
1222
1223
1224
1225
         else
1226
              //! if the load marker configuration button was clicked show a open file dialog
1227
              fileName = QFileDialog::getOpenFileName(nullptr, "Choose a previous saved marker configuration file
1228
      ", "", "marker configuratio files (*.xml);;All Files (*)");
1229
1230
              //! was cancel or abort clicked
1231
              if (fileName.length() == 0)
1232
              {
1233
                   //! if yes load the standard marker configuration
                  fileName = "markerStandard.xml";
1234
1235
1236
              //! open the selected marker configuration file
1237
1238
              FileStorage fs:
1239
              fs.open(fileName.toUtf8().constData(), FileStorage::READ);
1240
1241
              //! copy the values to the respective variables
1242
              fs["numberMarkers"] >> numberMarkers;
1243
1244
              //! inizialise vectors with correct length depending on the number of markers
1245
              list_points3d = std::vector<Point3d>(numberMarkers);
1246
              list_points2d = std::vector<Point2d>(numberMarkers);
1247
              list_points2dOld = std::vector<Point2d>(numberMarkers);
              list_points2dDifference = std::vector<double>(numberMarkers);
list_points2dProjected = std::vector<Point2d>(numberMarkers);
list_points2dUnsorted = std::vector<Point2d>(numberMarkers);
1248
1249
1250
1251
1252
              //! save the marker locations in the points3d vector
1253
              fs["list_points3d"] >> list_points3d;
1254
              fs.release();
1255
              commObj.addLog("Loaded marker configuration from file:");
1256
              commObj.addLog(fileName);
1257
1258
1259
1260
         //! Print out the number of markers and their position to the GUI
1261
         ss.str("");
         ss << "Number of Markers: " << numberMarkers << "\n";
12.62
         ss << "Marker 3D Points X, Y and Z [mm]: \n";
1263
1264
         for (int i = 0; i < numberMarkers; i++)</pre>
1265
      'ss << "Marker " << i + 1 << ":\t" << list_points3d[i].x << "\t" << list_points3d[i].y << "\t" << list_points3d[i].z << "\n";
1266
1267
         commObj.addLog(QString::fromStdString(ss.str()));
1268
1269
1270
         //! check if P3P algorithm can be enabled, it needs exactly 4 marker points to work
1271
          if (numberMarkers == 4)
1272
1273
              //! if P3P is possible, let the user choose which algorithm he wants but keep iterative active
              methodPNP = 0:
1274
1275
              commObj.enableP3P(true);
1276
1277
         else
1278
1279
              //! More (or less) marker than 4 loaded, P3P is not possible, hence user cant select P3P in GUI
1280
              methodPNP = 0:
              commObj.enableP3P(false);
1281
1282
              commObj.addLog("P3P algorithm disabled, only works with 4 markers.");
1283
1284
1285
          //! now display the marker configuration in the camera view
1286
         Mat cFrame(480, 640, CV_8UC3, Scalar(0, 0, 0));
1287
1288
          //! Set the camera pose parallel to the marker coordinate system
          Tvec.at<double>(0) = 0;
1289
1290
         Tvec.at < double > (1) = 0;
1291
          Tvec.at<double>(2) = 4500;
         Rvec.at<double>(0) = 0 * 3.141592653589 / 180.0;
Rvec.at<double>(1) = 0 * 3.141592653589 / 180.0;
1292
1293
         Rvec.at<double>(2) = -90. * 3.141592653589 / 180.0;
1294
1295
         projectPoints(list_points3d, Rvec, Tvec, cameraMatrix,
1296
      distCoeffs, list_points2dProjected);
1297
         for (int i = 0; i < numberMarkers; i++)</pre>
1298
              circle(cFrame, Point(list points2dProjected[i].x, list points2dProjected[i].y), 3, Scalar(255, 0, 0
1299
```

```
), 3);
1300 }
1301
1301
1302 projectCoordinateFrame(cFrame);
1303 QPixmap QPFrame;
1304 QPFrame = Mat2QPixmap(cFrame);
1305 commObj.changeImage(QPFrame);
1306 QCoreApplication::processEvents();
1307
1308 }
```

Here is the call graph for this function:



Here is the caller graph for this function:



3.1.2.11 main()

```
int main (
          int argc,
          char * argv[] )
```

main initialises the GUI and values for the marker position etc

First the GUI is set up with Signals and Slots, see Qt docu for how that works. Then some variables are initialized with arbitrary values. At last calibration and marker configuration etc. are loaded from xml files.


Parameters

| in | argc | is not used. |
|----|------|-------------------|
| in | argv | is also not used. |

Definition at line 156 of file main.cpp.

```
157 {
158
        QApplication a(argc, argv);
159
        RigidTrack w;
                     //!< show the GUI
160
         w.show();
161
         //! connect the Qt slots and signals for event handling
         QObject::connect(&commObj, SIGNAL(statusChanged(QString)), &w, SLOT(setStatus(QString)),
162
      Qt::DirectConnection);
163
        QObject::connect(&commObj, SIGNAL(imageChanged(QPixmap)), &w, SLOT(setImage(QPixmap)),
      Qt::DirectConnection);
        Object::connect(&commObj, SIGNAL(logAdded(OString)), &w, SLOT(setLog(OString)),
164
      Ot::DirectConnection);
165
        QObject::connect(&commObj, SIGNAL(logCleared()), &w, SLOT(clearLog(QString)),
      Qt::DirectConnection);
166
        QObject::connect(&commObj, SIGNAL(P3Penabled(bool)), &w, SLOT(enableP3P(bool)),
      Qt::DirectConnection);
        QObject::connect(&commObj, SIGNAL(progressUpdated(int)), &w, SLOT(progressUpdate(int)),
167
      Qt::DirectConnection);
168
169
         commObj.addLog("RigidTrack Version:");
170
         commObj.addLog(QString::number(_MSC_FULL_VER));
171
        commObj.addLog("Built on:");
         commObj.addLog(QString(__DATE__));
172
173
174
         //! initial guesses for position and rotation, important for Iterative Method!
175
         Tvec.at <double > (0) = 45;
176
         Tvec.at <double > (1) = 45;
         Tvec.at<double>(2) = 4500;
177
        Rvec.at<double>(1) = 0 * 3.141592653589 / 180.0;
Rvec.at<double>(1) = 0 * 3.141592653589 / 180.0;
178
179
        Rvec.at<double>(2) = -45 * 3.141592653589 / 180.0;
180
181
         //! Points that make up the marker CoSy axis system, hence one line in each axis direction
182
183
         coordinateFrame = std::vector<Point3d>(4);
         coordinateFrameProjected = std::vector<Point2d>(4);
184
        coordinateFrame[0] = cv::Point3d(0, 0, 0);
coordinateFrame[1] = cv::Point3d(300, 0, 0);
185
186
187
         coordinateFrame[2] = cv::Point3d(0, 300, 0);
188
         coordinateFrame[3] = cv::Point3d(0, 0, 300);
189
190
         position[0] = 1.1234;
                                      //! < \, {\rm set} \, \, {\rm position} \, \, {\rm initial} \, \, {\rm values}
                                      //!< set position initial values
        position[1] = 1.2345;
191
192
        position[2] = 1.3456;
                                      //!< set position initial values
193
194
         velocity[0] = 0.123;
                                    //!< set velocity initial values
195
         velocity[1] = 0.234;
                                    //!< set velocity initial values
        velocity[2] = 0.345;
196
                                    //! < \, {\rm set} \, \, {\rm velocity} \, \, {\rm initial} \, \, {\rm values}
197
        eulerAngles[0] = 1.002;
eulerAngles[1] = 1.003;
198
                                    //!< set initial euler angles to arbitrary values for testing
199
                                     //!< set initial euler angles to arbitrary values for testing
200
         eulerAngles[2] = 1.004;
                                    //!< set initial euler angles to arbitrary values for testing
201
202
         setHeadingOffset(0.0); //!< set the heading offset to 0</pre>
203
204
         ss.precision(4); //!< outputs in the log etc are limited to 3 decimal values
205
206
         loadCameraPosition(); //!< load the rotation matrix from camera CoSy to ground CoSy</pre>
        loadCalibration(0); //!< load the calibration file with the camera intrinsics
loadMarkerConfig(0); //!< load the standard marker configuration</pre>
207
208
209
         testAlgorithms(); //!< test the algorithms and their accuracy</pre>
210
211
         return a.exec();
212 }
```

Here is the call graph for this function:



3.1.2.12 Mat2QPixmap()

Convert an opency matrix that represents a picture to a Qt Pixmap object for the GUI.


Parameters

```
in src is the camera image represented as OpenCV matrix.
```

Definition at line 216 of file main.cpp.

```
217 {
218    QImage dest((const uchar *)src.data, src.cols, src.rows, src.step, QImage::Format_RGB888);
219    dest.bits(); //! enforce deep copy, see documentation
220    //! of QImage::QImage ( const uchar * data, int width, int height, Format format )
221    QPixmap pixmapDest = QPixmap::fromImage(dest);
222    return pixmapDest;
223 }
```

Here is the caller graph for this function:



3.1.2.13 projectCoordinateFrame()


Project the coordinate CoSy origin and axis direction of the marker CoSy with the rotation and translation of the object for visualization.

Parameters

| | in | pictureFrame | the image in which the CoSy frame should be pasted. | |
|--|----|--------------|---|--|
|--|----|--------------|---|--|

Definition at line 1081 of file main.cpp.

Here is the caller graph for this function:



3.1.2.14 sendDataUDP()

Send the position and attitude over UDP to every receiver, the safety receiver is handled on its own in the start

Tracking function because its send rate is less than 100 Hz.

Definition at line 1154 of file main.cpp.

```
1155 {
1156
         datagram.clear();
1157
         QDataStream out(&datagram, QIODevice::WriteOnly);
         out.setVersion(QDataStream::Qt_4_3);
1159
         out << (float)Position[0] << (float)Position[1] << (float)Position[2];</pre>
1160
         out << (float)Euler[0] << (float)Euler[1] << (float)Euler[2]; //! Roll Pitch Heading</pre>
1161
         udpSocketObject->writeDatagram(datagram,
      IPAdressObject, portObject);
1162
1163
         //! if second receiver is activated send it also the tracking data
         if (safety2Enable)
1166
             udpSocketSafety2->writeDatagram(datagram,
      IPAdressSafety2, portSafety2);
1167
1168
1169 }
```

Here is the caller graph for this function:



3.1.2.15 setHeadingOffset()

```
void setHeadingOffset ( double d )
```

Add a heading offset to the attitude for the case it is wanted by the user.

Parameters

```
in d denotes heading offset in degrees.
```

Definition at line 1122 of file main.cpp.

```
1123 {
          headingOffset = d;
1124
          d = d * 3.141592653589 / 180.0; //! Convert heading offset from degrees to rad
1125
1126
          //! Calculate rotation about x axis
1128
          Mat R_x = (Mat_<double>(3, 3) <<
               1, 0, 0,
1129
1130
               0, 1, 0,
1131
               0, 0, 1
1132
1133
1134
          //! Calculate rotation about y axis
         Mat R_y = (Mat_<double>(3, 3) <<
1, 0, 0,
0, 1, 0,
0, 0, 1
1135
1136
1137
1138
1139
               );
1140
1141
          //! Calculate rotation about z axis
         Mat R_z = (Mat_<double>(3, 3) <<
cos(d), -sin(d), 0,
sin(d), cos(d), 0,
0, 0, 1);
1142
1143
1144
1145
1146
1147
          //! Combined rotation matrix
1148
1149
          M_HeadingOffset = R_z * R_y * R_x;
1150 }
```

Here is the caller graph for this function:



3.1.2.16 setReference()

```
int setReference ( )
```

Determine the initial position of the object that serves as reference point or as ground frame origin. Computes the pose 200 times and then averages it. The position and attitude are from now on used as navigation CoSy.


Definition at line 595 of file main.cpp.

```
596 {
597
        //! initialize the variables with starting values
598
        gotOrder = false;
        posRef = 0;
599
600
        eulerRef = 0:
601
        RmatRef = 0;
        Rvec = RvecOriginal;
602
603
        Tvec = TvecOriginal;
604
605
        determineExposure();
606
        ss.str("");
607
608
        commObj.addLog("Started reference coordinate determination.");
609
610
        CameraLibrary_EnableDevelopment();
611
        //! Initialize Camera SDK
612
        CameraLibrary::CameraManager::X();
613
614
        //! At this point the Camera SDK is actively looking for all connected cameras and will initialize
615
        //! them on it's own.
616
617
        //! Get a connected camera ========
618
        CameraManager::X().WaitForInitialization();
619
        Camera *camera = CameraManager::X().GetCamera();
620
621
        //! If no device connected, pop a message box and exit ==--
622
        if (camera == 0)
623
62.4
            commObj.addLog("No camera found!");
625
            return 1;
626
627
628
        //! Determine camera resolution to size application window ==----
        int cameraWidth = camera->Width();
int cameraHeight = camera->Height();
629
630
631
        camera->GetDistortionModel(distModel);
        cv::Mat matFrame(cv::Size(cameraWidth, cameraHeight), CV_8UC1);
632
633
634
        //! Set camera mode to precision mode, it directly provides marker coordinates
635
        camera->SetVideoType(Core::PrecisionMode);
636
637
        //! Start camera output ==--
638
        camera->Start();
639
640
        //! Turn on some overlay text so it's clear things are
641
        //! working even if there is nothing in the camera's view. ===---
642
        //! Set some other parameters as well of the camera
643
        camera->SetTextOverlay(true);
644
        camera->SetFrameRate(intFrameRate);
645
        camera->SetIntensitv(intIntensity);
646
        camera->SetIRFilter(true);
647
        camera->SetContinuousIR(false);
648
        camera->SetHighPowerMode(false);
649
        //! sample some frames and calculate the position and attitude. then average those values and use that
650
       as zero position
651
        int numberSamples = 0;
        int numberToSample = 200;
652
        {\tt double\ projectionError\ =\ 0;\ //!<\ difference\ between\ the\ marker\ points\ as\ seen\ by\ the\ camera\ and\ the}
653
       projected marker points with Rvec and Tvec
654
655
        while (numberSamples < numberToSample)</pre>
656
657
             //! Fetch a new frame from the camera ===---
658
            Frame *frame = camera->GetFrame();
659
660
            if (frame)
661
                 //! Ok, we've received a new frame, lets do something
662
663
                 //! with it.
664
                 if (frame->ObjectCount() == numberMarkers)
665
                     //!for(int i=0; i<frame->ObjectCount(); i++)
for (int i = 0; i < numberMarkers; i++)</pre>
666
667
668
                         cObject *obj = frame->Object(i);
669
670
                          list_points2dUnsorted[i] = cv::Point2d(obj->X(), obj->Y());
671
672
673
                     if (gotOrder == false)
674
675
                         determineOrder();
676
677
678
                     //! sort the 2d points with the correct indices as found in the preceeding order
       determination algorithm
679
                     for (int w = 0; w < numberMarkers; w++)</pre>
```

```
680
                        list_points2d[w] = list_points2dUnsorted[
681
      pointOrderIndices[w]];
682
                    list_points2dOld = list_points2dUnsorted;
683
684
685
                    //!Compute the pose from the 3D-2D corresponses
686
                    solvePnP(list_points3d, list_points2d,
      cameraMatrix, distCoeffs, Rvec, Tvec, useGuess,
      methodPNP);
687
                    //! project the marker 3d points with the solution into the camera image CoSy and calculate
688
       difference to true camera image
689
                    projectPoints(list_points3d, Rvec, Tvec,
      cameraMatrix, distCoeffs, list_points2dProjected);
690
                    projectionError = norm(list_points2dProjected,
      list_points2d);
691
692
                    double maxValue = 0;
693
                    double minValue = 0;
694
                    minMaxLoc(Tvec.at < double > (2), & minValue, & maxValue);
695
696
                    if (maxValue > 10000 || minValue < 0)
697
                        ss.str("");
698
699
                        ss << "Negative z distance, thats not possible. Start the set zero routine again or
       restart Programm.";
700
                        commObj.addLog(QString::fromStdString(ss.str()));
701
                        frame->Release();
702
                        return 1:
703
                    }
704
705
                    if (projectionError > 3)
706
707
                        \verb|commObj.addLog("Reprojection error is bigger than 3 pixel. Correct marker| \\
       configuration loaded?\nMarker position measured precisely?");
708
                        frame->Release();
709
                        return 1;
710
711
                    if (norm(positionOld) - norm(Tvec) < 0.05) //!<Iterative Method needs time</pre>
712
       to converge to solution
713
714
                        add(posRef, Tvec, posRef);
                        add(eulerRef, Rvec, eulerRef); //!< That are not the values of yaw,
715
       roll and pitch yet! Rodriguez has to be called first.
716
                        numberSamples++;
                                            //!< one sample more :D
717
                        \verb|commObj.progressUpdate(numberSamples * 100 / numberToSample);|\\
718
719
                    positionOld = Tvec;
720
721
                    Mat cFrame(480, 640, CV_8UC3, Scalar(0, 0, 0));
722
                    for (int i = 0; i < numberMarkers; i++)</pre>
723
                        circle(cFrame, Point(list_points2d[i].x,
724
      list_points2d[i].y), 6, Scalar(0, 225, 0), 3);
725
726
                    projectCoordinateFrame(cFrame);
727
                    projectPoints(list_points3d, Rvec, Tvec,
      728
729
730
                        circle(cFrame, Point(list_points2d[i].x,
      list_points2d[i].y), 3, Scalar(225, 0, 0), 3);
731
732
                    drawPositionText(cFrame, position,
      eulerAngles, projectionError);
733
734
                    OPixmap OPFrame:
735
                    QPFrame = Mat2QPixmap(cFrame);
736
                     commObj.changeImage(QPFrame);
737
                    QCoreApplication::processEvents();
738
739
740
                frame->Release();
741
            }
742
743
        //! Release camera ==--
744
        camera->Release();
745
746
        \ensuremath{//!}\mbox{Divide} by the number of samples to get the mean of the reference position
747
        divide(posRef, numberToSample, posRef);
        divide(eulerRef, numberToSample, eulerRef); //!< eulerRef is here in Axis Angle
748
       notation
749
750
        Rodrigues(eulerRef, RmatRef);
                                                      //!< axis angle to rotation matrix
751
        //!-- Euler Angles, finally
```

```
getEulerAngles(RmatRef, eulerRef); //!< rotation matrix to euler</pre>
          ss.str("");
754
          ss << "RmatRef is:\n";
          ss << RmatRef << "\n";
755
          ss << "Reference Position is:\n";
ss << posRef << "[mm] \n";
ss << "Reference Euler Angles are:\n";
756
757
759
          ss << eulerRef << "[deg] \n";
760
        //! compute the difference between last obtained TVec and the average Value
//! When it is large the iterative method has not converged properly so it is advised to start the
setReference() function once again
double error = norm(posRef) - norm(Tvec);
761
762
763
764
         if (error > 5.0)
765
766
               ss << "Caution, distance between reference position and last position is: " << error << "\n Start
         the set zero routine once again.";
767
768
          commObj.addLog(QString::fromStdString(ss.str()));
769
          commObj.progressUpdate(0);
770
771 }
```

Here is the call graph for this function:



3.1.2.17 startTracking()

```
int startTracking ( )
```

Start the loop that fetches frames, computes the position etc and sends it to other computers. This function is the core of this program, hence the pose estimation is done here.

Definition at line 261 of file main.cpp.

```
261
                          {
262
263
        gotOrder = false; //! The order of points, hence which entry in list_points3d corresponds to
264
       which in list_points2d is not calculated yet
  Rvec = RvecOriginal; //! Use the value of Rvec that was set in main() as starting value
265
       for the solvePnP algorithm
         Tvec = TvecOriginal; //! Use the value of Tvec that was set in main() as starting value
266
       for the solvePnP algorithm
2.67
        GetLocalTime(&logDate); //! Get the current date and time to name the log file
268
         //! Concat the log file name as followed. The file is saved in the folder /logs in the Rigid Track
269
       installation folder
      logFileName = "./logs/positionLog_" + QString::number(logDate.wDay) + "_" +
QString::number(logDate.wMonth) + "_" + QString::number(logDate.wYear);
270
      logFileName += "_" + QString::number(logDate.wHour) + "_" + QString::number(logDate.wMinute) + "_" + QString::number(logDate.wSecond) + ".txt";
271
272
        logName = logFileName.toStdString(); //! Convert the QString to a standard string
273
274
        determineExposure(); //! Get the exposure where the right amount of markers is
275
2.76
         //! For OptiTrack Ethernet cameras, it's important to enable development mode if you
277
        //! want to stop execution for an extended time while debugging without disconnecting
278
        //! the Ethernet devices. Lets do that now:
279
280
        CameraLibrary_EnableDevelopment();
281
        CameraLibrary::CameraManager::X(); //! Initialize Camera SDK
282
283
        //! At this point the Camera SDK is actively looking for all connected cameras and will initialize
284
        //! them on it's own
285
286
         //! Get a connected camera
287
        CameraManager::X().WaitForInitialization();
288
        Camera *camera = CameraManager::X().GetCamera();
289
290
        //! If no camera can be found, inform user in message log and exit function
291
        if (camera == 0)
292
        {
293
             commObj.addLog("No camera found!");
294
             return 1;
295
296
297
         //! Determine camera resolution to size application window
298
         int cameraWidth = camera->Width();
299
         int cameraHeight = camera->Height();
300
        camera->SetVideoType(Core::PrecisionMode); //! Set the camera mode to precision mode, it used
301
       {\tt greyscale} \ {\tt imformation} \ {\tt for} \ {\tt marker} \ {\tt property} \ {\tt calculations}
302
303
        camera->Start(); //! Start camera output
304
305
         //! Turn on some overlay text so it's clear things are
306
         //! working even if there is nothing in the camera's view
307
        camera->SetTextOverlay(true);
308
        camera->SetExposure(intExposure);
                                                 //! Set the camera exposure
        camera->SetIntensity(intIntensity); //! Set the camera infrared LED intensity
309
        camera->SetFrameRate(intFrameRate); //! Set the camera framerate to 100 Hz
310
         camera->SetIRFilter(true); //! Enable the filter that blocks visible light and only passes infrared
311
       liaht
        camera->SetHighPowerMode(true); //! Enable high power mode of the LEDs
camera->SetContinuousIR(false); //! Disable continuous LED light
312
313
314
        camera->SetThreshold(intThreshold); //! Set threshold for marker detection
315
316
         //! Create a new matrix that stores the grayscale picture from the camera
317
        Mat matFrame = Mat::zeros(cv::Size(cameraWidth, cameraHeight), CV_8UC1);
318
        {\tt QPixmap} {\tt QPFrame;} //! {\tt QPixmap} is the corresponding {\tt Qt} class that saves images
         //! Matrix that stores the colored picture, hence marker points, coordinate frame and reprojected
319
       points
320
        Mat cFrame(480, 640, CV_8UC3, Scalar(0, 0, 0));
321
322
         int v = 0; //! Helper variable used to kick safety switch
323
         //! Variables for the min and max values that are needed for sanity checks
324
        double maxValue = 0;
325
        double minValue = 0;
326
         int framesDropped = 0; //! Ff a marker is not visible or accuracy is bad increase this counter
        double projectionError = 0; //! Equals the quality of the tracking
327
328
329
        setUpUDP(); //! Open sockets and ports for UDP communication
330
331
        if (safetyEnable) //! If the safety feature is enabled send the starting message
332
333
             //! Send enable message, hence send a 9 and then a 1
334
             data.setNum((int)(9));
335
             udpSocketSafety->write(data);
336
             data.setNum((int)(1));
337
             udpSocketSafety->write(data);
```


```
338
        }
339
340
        //! Fetch a new frame from the camera
        341
       subseeding frame so the time starts at {\tt 0} in the logs
342
        while (!gotTime) //! While no new frame is received loop
343
344
            Frame *frame = camera->GetFrame(); //! Get a new camera frame
345
            if (frame) //! There is actually a new frame
346
                timeFirstFrame = frame->TimeStamp(); //! Get the time stamp for the first frame.
347
       It is subtracted for the following frames
348
                                    //! Release the frame so the camera can continue
                frame->Release();
                gotTime = true; //! Exit the while loop
349
350
351
        }
352
        //! Now enter the main loop that processes each frame and computes the pose, sends it and logs stuff while (!exitRequested) //! Check if the user has not pressed "Stop Tracking" yet
353
354
355
356
357
            Frame *frame = camera->GetFrame(); //! Fetch a new frame from the camera
358
359
            if (frame) //! Did we got a new frame or does the camera still need more time
360
361
                framesDropped++; //! Increase by one, if everything is okay it is decreased at the end of the
       loop again
362
363
                //! Only use this frame it the right number of markers is found in the picture
364
                if (frame->ObjectCount() == numberMarkers)
365
366
                     //! Get the marker points in 2D in the camera image frame and store them in the
       list_points2dUnsorted vector
367
                    //! The order of points that come from the camera corresponds to the Y coordinate
368
                    for (int i = 0; i < numberMarkers; i++)</pre>
369
370
                        cObject *obj = frame->Object(i);
371
                        list_points2dUnsorted[i] = cv::Point2d(obj->X(), obj->Y());
372
373
374
                    if (gotOrder == false) //! Was the order already determined? This is false for the
       first frame and from then on true
375
376
                        determineOrder(); //! Now compute the order
377
378
379
                    //! Sort the 2d points with the correct indices as found in the preceeding order
       determination algorithm
380
                    for (int w = 0; w < numberMarkers; w++)</pre>
381
382
                        list_points2d[w] = list_points2dUnsorted[
      pointOrderIndices[w]]; //! pointOrderIndices was calculated in determineOrder()
383
384
                    list_points2dOld = list_points2dUnsorted;
385
386
                    //! The first time the 2D-3D corresspondence was determined with gotOrder was okay.
387
                    //! But this order can change as the object moves and the marker objects appear in a
388
                    //! different order in the frame->Object() array.
389
                     //! The solution is that: When a marker point (in the camera image, hence in 2D) was at
390
                     //! a position then it wont move that much from one frame to the other
391
                    //! So for the new frame we take a marker object and check which marker was closest this
       point
392
                    //! in the old image frame? This is probably the same (true) marker. And we do that for
       every other marker as well.
393
                    //! When tracking is good and no frames are dropped because of missing markers this should
       work every frame.
394
                    for (int j = 0; j < numberMarkers; j++)</pre>
395
                        minPointDistance = 5000; //! The sum of point distances is set to
396
       something unrealistic large
397
                        for (int k = 0; k < numberMarkers; k++)</pre>
398
399
                             //! Calculate N_2 norm of unsorted points minus old points
                            currentPointDistance = norm(
400
      list_points2dUnsorted[pointOrderIndices[j]] -
      list_points2dOld[k]);
401
                            //! If the norm is smaller than minPointDistance the correspondence is more likely
       to be correct
402
                            if (currentPointDistance <
      minPointDistance)
403
404
                                 //! Update the array that saves the new point order
405
                                minPointDistance =
      currentPointDistance;
406
                                 pointOrderIndicesNew[j] = k;
407
                            }
                        }
408
```

```
409
                    }
410
411
                    //! Now the new order is found, set the point order to the new value
412
                    for (int k = 0; k < numberMarkers; k++)</pre>
413
                         pointOrderIndices[k] = pointOrderIndicesNew[k];
414
                         list_points2d[k] = list_points2dUnsorted[
415
      pointOrderIndices[k]];
416
417
                     //! Save the unsorted position of the marker points for the next loop
418
                    list_points2dOld = list_points2dUnsorted;
419
420
421
                     //!Compute the object pose from the 3D-2D corresponses
422
                     solvePnP(list_points3d, list_points2d,
      cameraMatrix, distCoeffs, Rvec, Tvec, useGuess,
      methodPNP):
423
424
                    //! Project the marker 3d points with the solution into the camera image CoSy and calculate
       difference to true camera image
                    projectPoints(list_points3d, Rvec, Tvec,
425
      cameraMatrix, distCoeffs, list_points2dProjected);
                    projectionError = norm(list_points2dProjected,
426
      list_points2d); //! Difference of true pose and found pose
427
428
                     //! Increase the framesDropped variable if accuracy of tracking is too bad
429
                     if (projectionError > 5)
430
431
                         framesDropped++;
432
                     }
433
                    else
434
                    {
435
                         framesDropped = 0; //! Set number of subsequent frames dropped to zero because error
       is small enough and no marker was missing
436
437
438
                     //! Get the min and max values from TVec for sanity check
                    minMaxLoc(Tvec.at < double > (2), & minValue, & maxValue);
439
440
                    //! Sanity check of values. negative z means the marker CoSy is behind the camera, that's
441
       not possible.
442
                     if (minValue < 0)</pre>
443
444
                         commObj.addLog("Negative z distance, that is not possible. Start the set zero
       routine again or restart Program.");
445
                         frame->Release(); //! Release the frame so the camera can move on
446
                         camera->Release(); //! Release the camera
447
                        {\tt closeUDP} (); //! Close all UDP connections so the programm can be closed later
       on and no resources are locked
448
                        return 1; //! Exit the function
449
                     }
450
451
                     //! Next step is the transformation from camera CoSy to navigation CoSy
                    //! Compute the relative object position from the reference position to the current one //! given in the camera CoSy: f\ T_C^{NM} = Tvec - Tvec_{Ref} \ f
452
453
                    subtract(Tvec, posRef, position);
454
455
456
                     //! Transform the position from the camera CoSy to the navigation CoSy with INS alligned
       457
458
      position;
459
                    position = V;
                                     //! Position is the result of the preceeding calculation
                    position[2] *= invertZ; //! Invert Z if check box in GUI is activated,
460
       hence height above ground is considered
461
462
                     \ensuremath{//!} Realtive angle between reference orientation and current orientation
463
                    Rodrigues(Rvec, Rmat); //! Convert axis angle respresentation to ordinary rotation
       matrix
464
465
                     //! The difference of the reference rotation and the current rotation
466
                     //! \f$ R_{ NM } = M_{ NC } \times R_{ CM } \f$
467
                    Rmat = RmatRef.t() *Rmat;
468
                    //! Euler Angles, finally
getEulerAngles(Rmat, eulerAngles); //! Get the euler angles
469
470
       from the rotation matrix
471
                    eulerAngles[2] += headingOffset; //! Add the heading offset to the
       heading angle
472
                    //! Compute the velocity with finite differences. Only use is the log file. It is done here
473
       because the more precise time stamp can be used
                    frameTime = frame->TimeStamp() - timeOld; //! Time between the old frame
474
       and the current frame
475
                    timeOld = frame->TimeStamp();
                                                      //! Set the old frame time to the current one
                    velocity[0] = (position[0] - positionOld[0]) /
476
      frameTime; //! Calculate the x velocity with finite differences
```

```
477
                    velocity[1] = (position[1] - positionOld[1]) /
      frameTime; //! Calculate the y velocity with finite differences
                   velocity[2] = (position[2] - positionOld[2]) /
478
      479
      next frame velocity calcuation
480
481
                    //! Send position and Euler angles over WiFi with 100 Hz
482
                   sendDataUDP(position, eulerAngles);
483
484
                    //! Save the values in a log file, values are:
                    //! Time sinc tracking started Position Euler Angles
485
                                                                              Velocity
                    logfile.open(logName, std::ios::app); //! Open the log file, the folder is
486
       RigidTrackInstallationFolder/logs
487
                    logfile << frame->TimeStamp() - timeFirstFrame << ";" <</pre>
      488
                    logfile << velocity[0] << ";" << velocity[1] << ";" <<
489
      velocity[2] << "\n";</pre>
490
                    logfile.close(); //! Close the file to save values
491
492
               //! Check if the position and euler angles are below the allowed value, if yes send OKAY signal
493
       (1), if not send shutdown signal (0)
494
               //! Absolute x, y and z position in navigation CoSy must be smaller than the allowed distance
495
                if (safetyEnable)
496
497
                    if ((abs(position[0]) < safetyBoxLength && abs(position[1]) <</pre>
      safetyBoxLength && abs(position[2]) < safetyBoxLength))</pre>
498
499
                       //! Absolute Euler angles must be smaller than allowed value. Heading is not considered
500
                       if ((abs(eulerAngles[0]) < safetyAngle && abs(eulerAngles[1]) <</pre>
      safetyAngle))
501
                            \ensuremath{//!} Send the OKAY signal to the desired computer every 5th time
502
503
                           if (v == 5) {
504
                               data.setNum((int)(1));
505
                               udpSocketSafety->write(data); //! Send the 1
506
                               v = 0; //! reset the counter that is needed for decimation to every 5th time
      step
507
                           1
508
509
                       //! The euler angles of the object exceeded the allowed euler angles, send the shutdown
       signal (0)
510
                       else
511
                           data.setNum((int)(0)); //! Send the shutdown signal, a 0
512
513
                           udpSocketSafetv->write(data):
514
                           commObj.addLog("Object exceeded allowed Euler angles, shutdown signal sent."
      ); //! Inform the user
515
516
                       }
517
                   //! The position of the object exceeded the allowed position, shut the object down
518
519
                   else
520
                    {
521
                       data.setNum((int)(0)); //! Send the shutdown signal, a 0
522
                       udpSocketSafety->write(data);
                       commObj.addLog("Object left allowed area, shutdown signal sent."); //! Inform
523
      the user
524
525
                   }
526
               }
527
528
               //! Inform the user if tracking system is disturbed (marker lost or so) or error was too big
               if (framesDropped > 10)
529
530
531
                    if (safetyEnable) //! Also send the shutdown signal
532
533
                       data.setNum((int)(0)); //! Send the shutdown signal, a 0
534
                       udpSocketSafety->write(data);
535
536
                    commObj.addLog("Lost marker points or precision was bad!"); //! Inform the user
537
                    framesDropped = 0;
538
539
540
                //! Rasterize the frame so it can be shown in the GUI
               frame->Rasterize(cameraWidth, cameraHeight, matFrame.step,
541
      BACKBUFFER BITSPERPIXEL, matFrame.data);
542
543
                //! Convert the frame from greyscale as it comes from the camera to rgb color
544
               cvtColor(matFrame, cFrame, COLOR_GRAY2RGB);
545
546
               //! Project (draw) the marker CoSy origin into 2D and save it in the cFrame image
547
               projectCoordinateFrame(cFrame);
```

```
548
549
                 //! Project the marker points from 3D to the camera image frame (2d) with the computed pose
550
                 projectPoints(list_points3d, Rvec, Tvec,
      551
552
                     //! Draw a circle around the projected points so the result can be better compared to the
553
       real marker position
554
                    //! In the resulting picture those are the red dots
      circle(cFrame, Point(list_points2d[i].x,
list_points2d[i].y), 3, Scalar(225, 0, 0), 3);
555
556
               }
557
558
                 //! Write the current position, attitude and error values as text in the frame
559
                 drawPositionText(cFrame, position, eulerAngles, projectionError);
560
                 //! Send the new camera picture to the GUI and call the GUI processing routine
561
                 QPixmap QPFrame;
562
563
                 QPFrame = Mat2QPixmap(cFrame);
564
                 commObj.changeImage(QPFrame); //! Update the picture in the GUI
565
                 QCoreApplication::processEvents(); //! Give Qt time to handle everything
566
567
                 //! Release the camera frame to fetch the new one
568
                 frame->Release();
569
            }
570
        }
571
        //! User choose to stop the tracking, clean things up closeUDP(); //! Close the UDP connections so resources are deallocated camera->Release(); //! Release camera
572
573
574
575
        return 0:
576 }
```

Here is the call graph for this function:



Here is the caller graph for this function:



3.1.2.18 testAlgorithms()

void testAlgorithms () $\,$


Project some points from 3D to 2D and then check the accuracy of the algorithms. Mainly to generate something that can be shown in the camera view so the user knows everything loaded correctly.

Definition at line 952 of file main.cpp.

```
953 {
954
955
        int methodPNP;
956
957
        std::vector<Point2d> noise(numberMarkers);
958
       RvecOriginal = Rvec:
959
960
       TvecOriginal = Tvec;
961
962
        projectPoints(list_points3d, Rvec, Tvec, cameraMatrix,
      distCoeffs, list_points2dProjected);
963
964
        ss.str("");
       ss << "Unsorted Points 2D Projected \n";
ss << list_points2dProjected << "\n";
965
966
967
        commObj.addLog(QString::fromStdString(ss.str()));
968
969
        Mat cFrame(480, 640, CV_8UC3, Scalar(0, 0, 0));
970
        for (int i = 0; i < numberMarkers; i++)</pre>
971
            circle(cFrame, Point(list_points2dProjected[i].x, list_points2dProjected[i].y), 6, Scalar(0, 255, 0
972
      ), 3);
973
974
975
        projectCoordinateFrame(cFrame);
976
977
        ss.str("");
978
        ss << "=
979
        ss << "======Projected Points ======\n";
980
        ss << list_points2dProjected << "\n";
981
982
        randn(noise, 0, 0.5);
        add(list_points2dProjected, noise, list_points2dProjected);
983
984
        ss << "======= With Noise Points =======\n";
985
986
        ss << list_points2dProjected << "\n";</pre>
987
        commObj.addLog(QString::fromStdString(ss.str()));
988
989
990
        bool useGuess = true;
991
        _methodPNP = 0; //! < 0 = iterative 1 = EPNP 2 = P3P 4 = UPNP //! < not used
992
993
        solvePnP(list_points3d, list_points2dProjected, cameraMatrix,
      distCoeffs, Rvec, Tvec, useGuess, _methodPNP);
994
995
        ss.str("");
996
        ss << "====
997
        ss << "=======\n";
        ss << "rvec: " << "\n";
998
        ss << Rvec << "\n";
ss << "tvec: " << "\n";
999
1000
        ss << Tvec << "\n";
1001
1002
1003
         commObj.addLog(QString::fromStdString(ss.str()));
1004
1005
          _methodPNP = 1; //! < 0 = iterative 1 = EPNP 2 = P3P 4 = UPNP UPnP not used
1006
         Rvec = cv::Mat::zeros(3, 1, CV_64F);
         Tvec = cv::Mat::zeros(3, 1, CV_64F);
solvePnP(list_points3d, list_points2dProjected, cameraMatrix,
1007
1008
      distCoeffs, Rvec, Tvec, useGuess, _methodPNP);
1009
1010
         ss.str("");
1011
         =======\n":
         1012
         ss << "rvec: " << "\n";
1013
        ss << Rvec << "\n";
ss << "tvec: " << "\n";
1014
1015
1016
        ss << Tvec << "\n";
1017
     projectPoints(list_points3d, Rvec, Tvec, cameraMatrix,
distCoeffs, list_points2dProjected);
   for (int i = 0; i < numberMarkers; i++)</pre>
1018
1019
1020
1021
             circle(cFrame, Point(list_points2dProjected[i].x, list_points2dProjected[i].y), 3, Scalar(255, 0, 0
      ), 3);
1022
1023
         OPixmap OPFrame:
1024
         QPFrame = Mat2QPixmap(cFrame);
1025
         commObj.changeImage(QPFrame);
```

```
1026
         QCoreApplication::processEvents();
1027
         commObj.addLog(QString::fromStdString(ss.str()));
1028
         if (numberMarkers == 4)
1029
              methodPNP = 2; //!< 0 = iterative 1 = EPNP 2 = P3P 4 = UPNP //!< not used
1030
             Rvec = cv::Mat::zeros(3, 1, CV_64F);
Tvec = cv::Mat::zeros(3, 1, CV_64F);
1031
1032
1033
             solvePnP(list_points3d, list_points2dProjected,
      cameraMatrix, distCoeffs, Rvec, Tvec, useGuess, _methodPNP);
1034
             ss.str("");
1035
             ss << "----\n";
1036
             ss << "======\\n";
1037
1038
            ss << "rvec: " << "\n";
            ss << Rvec << "\n";
ss << "tvec: " << "\n";
1039
1040
             ss << Tvec << "\n";
1041
1042
1043
            projectPoints(list_points3d, Rvec, Tvec, cameraMatrix,
     distCoeffs, list_points2dProjected);
1044
            for (int i = 0; i < numberMarkers; i++)</pre>
1045
             {
1046
                 circle(cFrame, Point(list_points2dProjected[i].x, list_points2dProjected[i].y), 3, Scalar(255,
      0, 0), 3);
1047
1048
             double projectionError = norm(list_points2dProjected, list_points2d);
1049
             putText(cFrame, "Testing Algorithms Finished", cv::Point(5, 420), 1, 1, cv::Scalar(255, 255, 255));
1050
             drawPositionText(cFrame, position, eulerAngles, projectionError)
1051
1052
            QPixmap QPFrame;
1053
             QPFrame = Mat2QPixmap(cFrame);
1054
             commObj.changeImage(QPFrame);
1055
             QCoreApplication::processEvents();
1056
             commObj.addLog(QString::fromStdString(ss.str()));
1057
1058
1059
         _methodPNP = 4; //! < 0 = iterative 1 = EPNP 2 = P3P 4 = UPNP //! < not used
1060
         Rvec = cv::Mat::zeros(3, 1, CV_64F);
1061
        Tvec = cv::Mat::zeros(3, 1, CV_64F);
1062
         \verb|solvePnP| (list\_points3d, list\_points2dProjected, cameraMatrix, \\
     distCoeffs, Rvec, Tvec, useGuess, _methodPNP);
1063
1064
         ss.str("");
         ss << "==
1065
         ss << "==========
1066
                                       UPNP
                                               ======\n";
        ss << "rvec: " << "\n";
1067
        ss << "rvec: " << \n";
ss << Rvec << "\n";
ss << "tvec: " << "\n";
1068
1069
        ss << Tvec << "\n";
1070
1071
1072
        commObj.addLog(QString::fromStdString(ss.str()));
1073
        Rvec = RvecOriginal;
Tvec = TvecOriginal;
1074
1075
1076
1077 }
```

Here is the call graph for this function:



Here is the caller graph for this function:



3.1.3 Variable Documentation

3.1.3.1 commObj

```
commObject commObj
```

class that handles the communication from main.cpp to the GUI

Now declare variables that are used across the main.cpp file. Basically almost every variable used is declared here.

Definition at line 68 of file main.cpp.

3.1.3.2 Rmat

```
Mat Rmat = (cv::Mat_<double>(3, 1) << 0.0, 0.0, 0.0)
```

Rotation, translation etc. matrix for PnP results.

rotation matrix from camera CoSy to marker CoSy

Definition at line 95 of file main.cpp.


3.2 RigidTrack/main.h File Reference

Header file for main.cpp.

```
#include <fstream>
#include <windows.h>
#include <conio.h>
#include <tchar.h>
#include <stdio.h>
#include <iostream>
#include <stdarg.h>
#include <ctype.h>
#include <stdlib.h>
#include <gl/glu.h>
#include <sstream>
#include <thread>
#include <future>
#include <atomic>
#include "communication.h"
#include "RigidTrack.h"
#include <QtWidgets/QApplication>
#include <QUdpSocket>
#include "cameralibrary.h"
#include "modulevector.h"
#include "modulevectorprocessing.h"
#include "coremath.h"
#include <opencv\cv.h>
#include "opencv2\core.hpp"
#include "opencv2\calib3d.hpp"
#include <opencv2/imgproc/imgproc.hpp>
#include <opencv2/calib3d/calib3d.hpp>
#include <opencv2/highgui/highgui.hpp>
#include <opencv2\video\tracking.hpp>
Include dependency graph for main.h:
```



This graph shows which files directly or indirectly include this file:



Functions

- int startTracking ()
- · void startStopCamera ()

Start or stop the tracking depending on if the camera is currently running or not.

- int setReference ()
- int calibrateCamera ()

Start the camera calibration routine that computes the camera matrix and distortion coefficients.

- void loadCalibration (int method)
- · void testAlgorithms ()
- void projectCoordinateFrame (Mat pictureFrame)
- void setUpUDP ()

Open the UDP ports for communication.

- void setHeadingOffset (double d)
- void sendDataUDP (cv::Vec3d &Position, cv::Vec3d &Euler)
- void closeUDP ()
- void loadMarkerConfig (int method)
- void drawPositionText (cv::Mat &Picture, cv::Vec3d &Position, cv::Vec3d &Euler, double error)
- void loadCameraPosition ()
- int determineExposure ()
- void determineOrder ()
- int calibrateGround ()

Variables

int methodPNP

solvePNP algorithm $0 = iterative \ 1 = EPNP \ 2 = P3P \ 4 = UPNP //! < 4$ and 1 are the same and not implemented correctly by OpenCV

· bool safetyEnable

is the safety feature enabled

· bool safety2Enable

is the second receiver enabled

· double safetyBoxLength

length of the safety area cube in meters

· int safetyAngle

bank and pitch angle protection in degrees

• QHostAddress IPAdressObject

IPv4 adress of receiver 1.

· QHostAddress IPAdressSafety

IPv4 adress of safety receiver.

QHostAddress IPAdressSafety2

IPv4 adress of receiver 2.

int portObject

Port of receiver 1.

· int portSafety

Port of the safety receiver.

int portSafety2

Port of receiver 2.

int invertZ

dummy variable to invert Z direction on request

commObject commObj

class that handles the communication from main.cpp to the GUI

3.2.1 Detailed Description

Header file for main.cpp.

Author

Florian J.T. Wachter

Version

1.0

Date

April, 8th 2017

3.2.2 Function Documentation

3.2.2.1 calibrateGround()

```
int calibrateGround ( )
```

Get the pose of the camera w.r.t the ground calibration frame. This frame sets the navigation frame for later results. The pose is averaged over 200 samples and then saved in the file referenceData.xml. This routine is basically the same as setReference.


Definition at line 1563 of file main.cpp.

```
1564 {
         //! initialize the variables with starting values
1566
         gotOrder = false;
         posRef = 0;
1567
1568
         eulerRef = 0:
         RmatRef = 0;
1569
1570
         Rvec = RvecOriginal;
1571
         Tvec = TvecOriginal;
1572
1573
         determineExposure();
1574
1575
         ss.str("");
         commObj.addLog("Started ground calibration");
1576
1577
1578
         CameraLibrary_EnableDevelopment();
1579
         //! Initialize Camera SDK
1580
         CameraLibrary::CameraManager::X();
1581
1582
         //! At this point the Camera SDK is actively looking for all connected cameras and will initialize
1583
         //! them on it's own.
1584
1585
         //! Get a connected camera ======
1586
         CameraManager::X().WaitForInitialization();
1587
         Camera *camera = CameraManager::X().GetCamera();
1588
1589
         //! If no device connected, pop a message box and exit ==--
1590
         if (camera == 0)
1591
         {
1592
             commObj.addLog("No camera found!");
1593
             return 1;
1594
         }
1595
1596
         //! Determine camera resolution to size application window ==----
1597
         int cameraWidth = camera->Width();
1598
         int cameraHeight = camera->Height();
1599
         camera->GetDistortionModel(distModel);
1600
         cv::Mat matFrame(cv::Size(cameraWidth, cameraHeight), CV 8UC1);
1601
1602
         //! Set camera mode to precision mode, it directly provides marker coordinates
1603
         camera->SetVideoType(Core::PrecisionMode);
1604
1605
         //! Start camera output ==--
1606
         camera->Start();
1607
1608
           '! Turn on some overlay text so it's clear things are
1609
         //! working even if there is nothing in the camera's view. ===---
1610
         //! Set some other parameters as well of the camera
1611
         camera->SetTextOverlay(true);
1612
         camera->SetFrameRate(intFrameRate);
1613
         camera->SetIntensity(intIntensity);
1614
         camera->SetIRFilter(true);
1615
         camera->SetContinuousIR(false);
1616
         camera->SetHighPowerMode(false);
1617
1618
         //! sample some frames and calculate the position and attitude. then average those values and use that
       as zero position
1619
         int numberSamples = 0;
         int numberToSample = 200;
1620
1621
         double projectionError = 0;
1622
1623
         while (numberSamples < numberToSample)</pre>
1624
             //! Fetch a new frame from the camera ===---
1625
1626
             Frame *frame = camera->GetFrame();
1627
1628
                (frame)
1629
                  //! Ok, we've received a new frame, lets do something
1630
1631
                 //! with it.
1632
                  if (frame->ObjectCount() == numberMarkers)
1633
```

```
1634
                      //!for(int i=0; i<frame->ObjectCount(); i++)
1635
                      for (int i = 0; i < numberMarkers; i++)</pre>
1636
1637
                          cObject *obj = frame->Object(i);
1638
                          list_points2dUnsorted[i] = cv::Point2d(obj->X(), obj->Y());
1639
1640
1641
                      if (gotOrder == false)
1642
1643
                          determineOrder();
1644
1645
                      //! sort the 2d points with the correct indices as found in the preceeding order
1646
       determination algorithm
1647
                      for (int w = 0; w < numberMarkers; w++)</pre>
1648
1649
                          list_points2d[w] = list_points2dUnsorted[
      pointOrderIndices[w]];
1650
1651
                      list_points2dOld = list_points2dUnsorted;
1652
1653
                      //!Compute the pose from the 3D-2D corresponses
1654
                      solvePnP(list_points3d, list_points2d,
      cameraMatrix, distCoeffs, Rvec, Tvec, useGuess,
      methodPNP);
1655
                      //! project the marker 3d points with the solution into the camera image CoSy and calculate
1656
       difference to true camera image
1657
                      projectPoints(list_points3d, Rvec, Tvec,
      cameraMatrix, distCoeffs, list_points2dProjected);
1658
                     projectionError = norm(list_points2dProjected,
      list_points2d);
1659
1660
                      if (projectionError > 3)
1661
                          commObj.addLog("Reprojection error is bigger than 3 pixel. Correct marker
1662
       configuration loaded?\nMarker position measured precisely?");
1663
                         frame->Release();
1664
                          return 1;
1665
1666
1667
                      double maxValue = 0;
                      double minValue = 0;
1668
                      minMaxLoc(Tvec.at < double > (2), & minValue, & maxValue);
1669
1670
1671
                      if (maxValue > 10000 || minValue < 0)</pre>
1672
1673
1674
                          commObj.addLog("Negative z distance, thats not possible. Start the set zero
1675
       routine again and check marker configurations.");
1676
                          frame->Release();
1677
                          return 1;
1678
                      }
1679
                      if (norm(positionOld) - norm(Tvec) < 0.05) //!<Iterative Method needs time</pre>
1680
       to converge to solution
1681
                     {
1682
                          add(posRef, Tvec, posRef);
                          add(eulerRef, Rvec, eulerRef); //!< That are not the values of yaw,
1683
       roll and pitch yet! Rodriguez has to be called first. 
 numberSamples++; \qquad //! <-- \ one \ sample \ more :D
1684
1685
                          commObj.progressUpdate(numberSamples * 100 / numberToSample);
1686
                     positionOld = Tvec;
1687
1688
1689
                      Mat cFrame(480, 640, CV_8UC3, Scalar(0, 0, 0));
1690
                      for (int i = 0; i < numberMarkers; i++)</pre>
1691
1692
                          circle(cFrame, Point(list_points2d[i].x,
      list_points2d[i].y), 6, Scalar(0, 225, 0), 3);
1693
                      projectCoordinateFrame(cFrame);
1694
      1695
1696
1697
                      {
1698
                         circle(cFrame, Point(list_points2d[i].x,
      list_points2d[i].y), 3, Scalar(225, 0, 0), 3);
1699
1700
1701
                      QPixmap QPFrame;
                      QPFrame = Mat2QPixmap(cFrame);
1702
1703
                      commObj.changeImage(QPFrame);
1704
                      QCoreApplication::processEvents();
1705
1706
                 }
```

```
frame->Release();
1708
1709
                                //! Release camera ==--
1710
1711
                                camera->Release();
1712
1713
                                 \label{eq:continuous} \parbox{0.05\line Model}{\parbox{0.05\line Mode
1714
                                divide(posRef, numberToSample, posRef);
1715
                                divide(eulerRef, numberToSample, eulerRef); //!< eulerRef is here in Axis Angle
1716
1717
                                Rodrigues(eulerRef, RmatRef);
                                                                                                                                                                                             //!< axis angle to rotation matrix
1718
1719
                                getEulerAngles(RmatRef, eulerRef); //!< rotation matrix to euler</pre>
                               ss.str("");
ss << "RmatRef is:\n";
ss << RmatRef << "\n";</pre>
1720
1721
1722
                              ss << "Reference Position is:\n";
ss << posRef << "[mm] \n";
1723
1724
1725
                                ss << "Reference Euler angles are:\n";
1726
                                ss << eulerRef << "[deg] \n";
1727
                   //! Save the obtained calibration coefficients in a file for later use
    QString fileName = QFileDialog::getSaveFileName(nullptr, "Save ground calibration file", "
referenceData.xml", "Calibration File (*.xml);;All Files (*)");
    FileStorage fs(fileName.toUtf8().constData(), FileStorage::WRITE);
    fs << "M_NC" << RmatRef;
    fs << "eulerRef;
    fs << "eulerRef" << eulerRef;</pre>
1728
1729
1730
1731
1732
                                strBuf = fs.releaseAndGetString();
1733
                                commObj.changeStatus(QString::fromStdString(strBuf));
1734
                                commObj.addLog("Saved ground calibration!");
1735
1736
                                commObj.progressUpdate(0);
1737
                                return 0;
1738 }
```

Here is the call graph for this function:



3.2.2.2 closeUDP()


void closeUDP ()

Close the UDP ports again to release network interfaces etc. If this is not done the network resources are still occupied and the program can't exit properly.

Definition at line 1173 of file main.cpp.

```
1174 {
1175
         //! check if the socket is open and if yes close it
1176
         if (udpSocketObject->isOpen())
1177
1178
             udpSocketObject->close();
1179
1180
1181
        if (udpSocketSafety->isOpen())
1182
1183
             udpSocketSafety->close();
1184
1185
         if (udpSocketSafety2->isOpen())
1186
1187
1188
             udpSocketSafety2->close();
1189
1190
         commObj.addLog("Closed all UDP ports.");
1191 }
```

Here is the caller graph for this function:



3.2.2.3 determineExposure()

```
int determineExposure ( )
```

Get the optimal exposure for the camera. For that find the minimum and maximum exposure were the right number of markers are detected. Then the mean of those two values is used as exposure.


Definition at line 1362 of file main.cpp.

```
1363 {
1364
         //! For OptiTrack Ethernet cameras, it's important to enable development mode if you
1365
         //! want to stop execution for an extended time while debugging without disconnecting
         //! the Ethernet devices. Lets do that now:
1366
1367
1368
         CameraLibrary_EnableDevelopment();
1369
1370
         //! Initialize Camera SDK ==--
1371
         CameraLibrary::CameraManager::X();
1372
1373
         //! At this point the Camera SDK is actively looking for all connected cameras and will initialize
1374
         //! them on it's own.
1375
1376
         //! Get a connected camera ==
1377
         {\tt CameraManager::X().WaitForInitialization();}
1378
         Camera *camera = CameraManager::X().GetCamera();
1379
1380
         //! If no device connected, pop a message box and exit ==--
1381
         if (camera == 0)
```

```
1382
               {
1383
                      commObj.addLog("No camera found!");
1384
                      return 1;
1385
               }
1386
               //! Determine camera resolution to size application window ==----
1387
1388
               int cameraWidth = camera->Width();
1389
               int cameraHeight = camera->Height();
1390
               camera->SetVideoType(Core::PrecisionMode); //! set the camera mode to precision mode, it used
1391
           greyscale imformation for marker property calculations
1392
1393
                                                                                          //! Start camera output ==--
1394
               camera->Start();
1395
1396
               //! Turn on some overlay text so it's clear things are
1397
               //! working even if there is nothing in the camera's view. ===---
1398
               camera->SetTextOverlay(true);
1399
               camera->SetExposure(intExposure);
                                                                               //! set the camera exposure
               camera->SetIntensity(intIntensity); //! set the camera infrared LED intensity
camera->SetFrameRate(intFrameRate); //! set the camera framerate to 100 Hz
1400
1401
1402
               camera->SetIRFilter(true); //! enable the filter that blocks visible light and only passes infrared
            light
               camera->SetHighPowerMode(true); //! enable high power mode of the leds
camera->SetContinuousIR(false); //! enable continuous LED light
1403
1404
               camera->SetThreshold(intThreshold); //! set threshold for marker detection
1405
1406
1407
               //!set exposure such that num markers are visible
1408
               int numberObjects = 0; //! Number of objects (markers) found in the current picture with the given
           exposure
1409
               int minExposure = 1; //! exposure when objects detected the first time is numberMarkers int maxExposure = 480; //! exposure when objects detected is first time numberMarkers+1
1410
1411
               intExposure = minExposure; //! set the exposure to the smallest value possible
               int numberTries = 0; //! if the markers arent found after numberTries then there might be no markers
1412
            at all in the real world
1413
                                                        //! Determine minimum exposure, hence when are numberMarkers objects detected
1414
               camera->SetExposure(intExposure);
1415
1416
               while (numberObjects != numberMarkers && numberTries < 48)</pre>
1417
1418
                       //! get a new camera frame
                      Frame *frame = camera->GetFrame();
if (frame) //! frame received
1419
1420
1421
                             numberObjects = frame->ObjectCount();  //! how many objects are detected in the image
1422
1423
                             if (numberObjects == numberMarkers) { minExposure =
          intExposure; frame->Release(); break; } //! if the right amount if markers is found, exit while
1424
                             //! not the right amount of markers was found so increase the exposure and try again
1425
                             numberTries++;
                             intExposure += 10;
1426
1427
                             camera->SetExposure(intExposure);
                             ss.str("");
1428
                            ss << "Exposure: " << intExposure << "\t";
ss << "Objects found: " << numberObjects;
1429
1430
                             commObj.addLog(QString::fromStdString(ss.str()));
1431
1432
                             frame->Release();
1433
                      }
1434
1435
               //! Now determine maximum exposure, hence when are numberMarkers+1 objects detected
1436
               numberTries = 0; //! if the markers arent found after numberTries then there might be no markers at
1437
           all in the real world
1438
              intExposure = maxExposure;
1439
               camera->SetExposure(intExposure);
1440
               numberObjects = 0;
1441
               while (numberObjects != numberMarkers && numberTries < 48)</pre>
1442
1443
                      Frame *frame = camera->GetFrame();
1444
                      if (frame)
1445
1446
                             \verb|numberObjects| = frame->ObjectCount(); //! | how many objects | are detected in the image | leaves | frame->ObjectCount(); | frame->ObjectCount();
          if (numberObjects == numberMarkers) { maxExposure =
intExposure; frame->Release(); break; } //! if the right amount if markers is found, exit while
1447
           100p
1448
1449
                             //! not the right amount of markers was found so decrease the exposure and try again
1450
                             intExposure -= 10;
1451
                             numberTries++:
                             camera->SetExposure(intExposure);
1452
                             ss.str("");
1453
                             ss << "Exposure: " << intExposure << "\t";
1454
                             ss << "Objects found: " << numberObjects;
1455
1456
                             commObj.addLog(QString::fromStdString(ss.str()));
1457
                             frame->Release();
1458
                      }
              }
1459
```

```
1460
1461
          //! set the exposure to the mean of min and max exposure determined
1462
          camera->SetExposure((minExposure + maxExposure) / 2.0);
1463
1464
          //! and now check if the correct amount of markers is detected with that new value
1465
          while (1)
1466
          {
1467
               Frame *frame = camera->GetFrame();
1468
                if (frame)
1469
                    numberObjects = frame->ObjectCount(); //! how many objects are detected in the image
if (numberObjects != numberMarkers) //! are all markers and not more or less
1470
1471
        detected in the image
1472
1473
                         frame->Release();
1474
                         \verb|commObj.addLog("Was not able to detect the right amount of markers.");\\
1475
                         //! Release camera ==-
1476
                         camera->Release();
1477
                         return 1;
1478
1479
                    else //! all markers and not more or less are found
1480
                         frame->Release();
1481
                         intExposure = (minExposure + maxExposure) / 2.0;
commObj.addLog("Found the correct number of markers.");
commObj.addLog("Exposure set to:");
1482
1483
1484
1485
                         commObj.addLog(QString::number(intExposure));
1486
1487
1488
               }
1489
          }
1490
1491
          camera->Release();
1492
          return 0;
1493
1494 }
```

Here is the caller graph for this function:



3.2.2.4 determineOrder()


```
void determineOrder ( )
```

Compute the order of the marker points in 2D so they are the same as in the 3D array. Hence marker 1 must be in first place for both, list_points2d and list_points3d.

Definition at line 1498 of file main.cpp.

```
1499 {
1500
                           //! determine the 3D-2D correspondences that are crucial for the PnP algorithm
1501
                          //! Try every possible correspondence and solve PnP
                          //! Then project the 3D marker points into the 2D camera image and check the difference
1502
1503
                          //! between projected points and points as seen by the camera
                          //! the corresponce with the smallest difference is probably the correct one
1504
1505
1506
                                       //! the difference between true 2D points and projected points is super big
1507
                          minPointDistance = 5000;
1508
                          std::sort(pointOrderIndices, pointOrderIndices + 4);
1509
1510
                          //! now try every possible permutation of correspondence
1511
                                       //! reset the starting values for solvePnP
1512
1513
                                       Rvec = RvecOriginal;
                                      Tvec = TvecOriginal;
1514
1515
                                      //! sort the 2d points with the current permutation
1516
                                      for (int m = 0; m < numberMarkers; m++)</pre>
1517
1518
                                      {
                                                  list_points2d[m] = list_points2dUnsorted[
1519
                pointOrderIndices[m]];
1520
                                   }
1521
1522
                                      //! Call solve PNP with P3P since its more robust and sufficient for start value determination
1523
                                       solvePnP(list_points3d, list_points2d,
                 cameraMatrix, distCoeffs, Rvec, Tvec, useGuess, SOLVEPNP_P3P);
1524
1525
                                       //! set the current difference of all point correspondences to zero
1526
                                      currentPointDistance = 0;
1527
1528
                                      //! project the 3D points with the solvePnP solution onto 2D
                                      projectPoints(list_points3d, Rvec, Tvec,
1529
                cameraMatrix, distCoeffs, list_points2dProjected);
1530
                                       //! now compute the absolute difference (error)
1531
1532
                                      for (int n = 0; n < numberMarkers; n++)</pre>
1533
1534
                                                  currentPointDistance += norm(list_points2d[n] -
                 list_points2dProjected[n]);
1535
                                     }
1536
                                      //! if the difference with the current permutation is smaller than the smallest value till now
1537
1538
                                      //! it is probably the more correct permutation
                                       if (currentPointDistance < minPointDistance)</pre>
1539
1540
1541
                                                 minPointDistance = currentPointDistance;
                                                                                                                                                                             //!< set the
                    smallest value of difference to the current one % \left( 1\right) =\left( 1\right) \left( 1\right
1542
                                                 for (int b = 0; b < numberMarkers; b++)</pre>
                                                                                                                                                                             //!< now safe the better permutation
1543
                                                  {
1544
                                                             pointOrderIndicesNew[b] = pointOrderIndices[b];
1545
1546
                                      }
1547
1548
1549
 1550
                          //! try every permutation
1551
                           while (std::next_permutation(pointOrderIndices,
                 pointOrderIndices + 4));
1552
1553
                           //! now that the correct order is found assign it to the indices array
1554
                          for (int w = 0; w < numberMarkers; w++)</pre>
1555
1556
                                      pointOrderIndices[w] = pointOrderIndicesNew[w];
1557
                          gotOrder = true;
1558
1559 }
```

Here is the caller graph for this function:



3.2.2.5 drawPositionText()

Draw the position, attitude and reprojection error in the picture.


Parameters

| in | Picture | is the camera image in OpenCV matrix format. | |
|----|----------|--|--|
| in | Position | is the position of the tracked object in navigation CoSy. | |
| in | Euler | are the Euler angles with respect to the navigation frame. | |
| in | error | is the reprojection error of the pose estimation. | |

Definition at line 1315 of file main.cpp.

```
1316 {
         ss.str("");
ss << "X: " << Position[0] << " m";</pre>
1317
1318
1319
         putText(Picture, ss.str(), cv::Point(200, 440), 1, 1, cv::Scalar(255, 255, 255));
1320
         ss.str("");
ss << "Y: " << Position[1] << " m";</pre>
1321
1322
         putText(Picture, ss.str(), cv::Point(200, 455), 1, 1, cv::Scalar(255, 255, 255));
1323
1324
         ss.str("");
ss << "Z: " << Position[2] << " m";</pre>
1325
1326
1327
         putText(Picture, ss.str(), cv::Point(200, 470), 1, 1, cv::Scalar(255, 255, 255));
1328
         ss.str("");
1329
1330
         ss << "Heading: " << Euler[2] << " deg";
1331
         putText(Picture, ss.str(), cv::Point(350, 440), 1, 1, cv::Scalar(255, 255, 255));
1332
         ss.str("");
ss << "Pitch: " << Euler[1] << " deg";</pre>
1333
1334
1335
         putText(Picture, ss.str(), cv::Point(350, 455), 1, 1, cv::Scalar(255, 255, 255));
1336
1337
          ss.str("");
```

Here is the caller graph for this function:



3.2.2.6 loadCalibration()

Load a previously saved camera calibration from a file.

Parameters

| in | method | whether or not load the camera calibration from calibration.xml. If ==0 then yes, if != 0 then | |
|----|--------|--|--|
| | | the user select a different file. | |

Definition at line 923 of file main.cpp.

```
923
924
925
        QString fileName;
926
         if (method == 0)
927
             fileName = "calibration.xml";
928
929
930
        else
931
932
              fileName = QFileDialog::getOpenFileName(nullptr, "Choose a previous saved calibration file", "", "
      Calibration Files (*.xml);;All Files (*)");
933
             if (fileName.length() == 0)
934
             {
935
                  fileName = "calibration.xml";
             }
936
937
        FileStorage fs;
938
939
        fs.open(fileName.toUtf8().constData(), FileStorage::READ);
        fs.Topen(fileName.tootib().constData
fs["CameraMatrix"] >> cameraMatrix;
fs["DistCoeff"] >> distCoeffs;
940
941
942
        commObj.addLog("Loaded calibration from file:");
```

```
943 commObj.addLog(fileName);

944 ss.str("");

945 ss << "\nCamera Matrix is" << "\n" << cameraMatrix << "\n";

946 ss << "\nDistortion Coefficients are" << "\n" << distCoeffs << "\n";

947 commObj.addLog(QString::fromStdString(ss.str()));

948 }
```

Here is the caller graph for this function:



3.2.2.7 loadCameraPosition()

```
void loadCameraPosition ( )
```

Load the rotation matrix from camera CoSy to ground CoSy It is determined during calibrateGround() and stays the same once the camera is mounted and fixed.

Definition at line 1348 of file main.cpp.

Here is the caller graph for this function:



3.2.2.8 loadMarkerConfig()

Load a marker configuration from file. This file has to be created by hand, use the standard marker configuration file as template.

Parameters

in *method*


whether or not load the configuration from the markerStandard.xml. If ==0 load it, if != 0 let the user select a different file.

Definition at line 1195 of file main.cpp.

```
1196 {
1197
         QString fileName;
         //! during start up of the programm load the standard marker configuration
1198
1199
         if (method == 0)
1200
1201
              //! open the standard marker configuration file
1202
              FileStorage fs;
1203
              fs.open("markerStandard.xml", FileStorage::READ);
1204
1205
              //! copy the values to the respective variables
1206
              fs["numberMarkers"] >> numberMarkers;
1207
1208
              //! inizialise vectors with correct length depending on the number of markers
1209
              list_points3d = std::vector<Point3d>(numberMarkers);
              list_points2d = std::vector<Point2d>(numberMarkers);
1210
1211
              list_points2dOld = std::vector<Point2d>(numberMarkers);
1212
              list_points2dDifference = std::vector<double>(
      numberMarkers);
1213
              list_points2dProjected = std::vector<Point2d>(
      numberMarkers);
              list_points2dUnsorted = std::vector<Point2d>(
1214
      numberMarkers);
1215
1216
              //! save the marker locations in the points3d vector
1217
              fs["list_points3d"] >> list_points3d;
1218
              fs.release();
              commObj.addLog("Loaded marker configuration from file:");
1219
              commObj.addLog(fileName);
1220
1221
1222
1223
1224
1225
         else
1226
              //! if the load marker configuration button was clicked show a open file dialog
1227
              fileName = QFileDialog::getOpenFileName(nullptr, "Choose a previous saved marker configuration file
1228
      ", "", "marker configuratio files (*.xml);; All Files (*)");
1229
1230
              //! was cancel or abort clicked
1231
              if (fileName.length() == 0)
1232
              {
                  //! if yes load the standard marker configuration
fileName = "markerStandard.xml";
1233
1234
1235
1236
1237
              //! open the selected marker configuration file
1238
              FileStorage fs;
1239
              fs.open(fileName.toUtf8().constData(), FileStorage::READ);
1240
              //! copy the values to the respective variables
1241
1242
              fs["numberMarkers"] >> numberMarkers;
1243
1244
              //! inizialise vectors with correct length depending on the number of markers
1245
              list_points3d = std::vector<Point3d>(numberMarkers);
              list_points2d = std::vector<Point2d>(numberMarkers);
1246
1247
              list_points2dOld = std::vector<Point2d>(numberMarkers);
              list_points2dDifference = std::vector<double>(numberMarkers);
list_points2dProjected = std::vector<Point2d>(numberMarkers);
list_points2dUnsorted = std::vector<Point2d>(numberMarkers);
1248
1249
1250
1251
1252
              //! save the marker locations in the points3d vector
1253
              fs["list_points3d"] >> list_points3d;
1254
              fs.release();
              commObj.addLog("Loaded marker configuration from file:");
1255
1256
              commObj.addLog(fileName);
1257
1258
         }
1259
1260
          //! Print out the number of markers and their position to the GUI
1261
         ss.str("");
         ss << "Number of Markers: " << numberMarkers << "\n";
1262
         ss << "Marker 3D Points X, Y and Z [mm]: \n";
1263
1264
         for (int i = 0; i < numberMarkers; i++)</pre>
1265
```

```
ss << "Marker" << i + 1 << ":\t" << list_points3d[i].x << "\t" << list_points3d[i].y << "\t" << li>x <= "\t" << li>x
1266
              list_points3d[i].z << "\n";
1267
1268
                    commObj.addLog(QString::fromStdString(ss.str()));
1269
                    //! check if P3P algorithm can be enabled, it needs exactly 4 marker points to work
1270
1271
                    if (numberMarkers == 4)
1272
1273
                              //! if P3P is possible, let the user choose which algorithm he wants but keep iterative active
1274
                              methodPNP = 0;
                              commObj.enableP3P(true);
1275
1276
1277
                    else
1278
1279
                              //! More (or less) marker than 4 loaded, P3P is not possible, hence user cant select P3P in GUI
1280
                              methodPNP = 0;
                              commObj.enableP3P(false);
1281
                              commObj.addLog("P3P algorithm disabled, only works with 4 markers.");
1282
1283
1284
1285
                     //! now display the marker configuration in the camera view
1286
                    Mat cFrame(480, 640, CV_8UC3, Scalar(0, 0, 0));
1287
                    \ensuremath{//!} Set the camera pose parallel to the marker coordinate system
1288
1289
                    Tvec.at < double > (0) = 0;
1290
                    Tvec.at < double > (1) = 0;
                    Tvec.at<double>(2) = 4500;
1291
                    Rvec.at<double>(0) = 0 * 3.141592653589 / 180.0;
Rvec.at<double>(1) = 0 * 3.141592653589 / 180.0;
1292
1293
1294
                    Rvec.at<double>(2) = -90. * 3.141592653589 / 180.0;
1295
1296
                    projectPoints(list_points3d, Rvec, Tvec, cameraMatrix,
             distCoeffs, list_points2dProjected);
1297
                     for (int i = 0; i < numberMarkers; i++)</pre>
1298
                              circle(cFrame, Point(list_points2dProjected[i].x, list_points2dProjected[i].y), 3, Scalar(255, 0, 0
1299
             ), 3);
1300
1301
1302
                    projectCoordinateFrame(cFrame);
1303
                    QPixmap QPFrame;
                    QPFrame = Mat2QPixmap(cFrame);
1304
1305
                    commObj.changeImage(OPFrame);
1306
                    QCoreApplication::processEvents();
1307
1308 }
```

Here is the call graph for this function:



Here is the caller graph for this function:




3.2.2.9 projectCoordinateFrame()

Project the coordinate CoSy origin and axis direction of the marker CoSy with the rotation and translation of the object for visualization.

Parameters

Definition at line 1081 of file main.cpp.

Here is the caller graph for this function:



3.2.2.10 sendDataUDP()

Send the position and attitude over UDP to every receiver, the safety receiver is handled on its own in the start

Tracking function because its send rate is less than 100 Hz.

Definition at line 1154 of file main.cpp.

```
1155 {
1156
          datagram.clear();
1157
          QDataStream out(&datagram, QIODevice::WriteOnly);
1158
          out.setVersion(QDataStream::Qt_43);
          out << (float)Position[0] << (float)Position[1] << (float)Position[2];
out << (float)Euler[0] << (float)Euler[1] << (float)Euler[2]; //! Roll Pitch Heading</pre>
1159
1160
          udpSocketObject->writeDatagram(datagram,
1161
      IPAdressObject, portObject);
1162
1163
          //! if second receiver is activated send it also the tracking data
1164
          if (safety2Enable)
1165
         {
              udpSocketSafety2->writeDatagram(datagram,
1166
      IPAdressSafety2, portSafety2);
1167
1168
1169 }
```

Here is the caller graph for this function:



3.2.2.11 setHeadingOffset()

```
void setHeadingOffset ( double d )
```

Add a heading offset to the attitude for the case it is wanted by the user.

Parameters

| in d denotes heading offset in degree | s. |
|---------------------------------------|----|
|---------------------------------------|----|

Definition at line 1122 of file main.cpp.

```
1123 {
1124
         headingOffset = d;
         d = d \star 3.141592653589 / 180.0; //! Convert heading offset from degrees to rad
1125
1126
         //! Calculate rotation about x axis
1127
1128
        1, 0, 0,
0, 1, 0,
1129
1130
1131
             0, 0, 1
1132
1133
         //! Calculate rotation about y axis
1134
1135
         Mat R_y = (Mat_{double})(3, 3) <<
             1, 0, 0,
0, 1, 0,
1136
1137
1138
             0, 0, 1
1139
             );
1140
1141
         //! Calculate rotation about z axis
1142
        Mat R_z = (Mat_{double})(3, 3) <<
```

Here is the caller graph for this function:



3.2.2.12 setReference()

```
int setReference ( )
```

Determine the initial position of the object that serves as reference point or as ground frame origin. Computes the pose 200 times and then averages it. The position and attitude are from now on used as navigation CoSy.

Definition at line 595 of file main.cpp.


```
596 {
597
        //! initialize the variables with starting values
598
        gotOrder = false;
        posRef = 0;
599
600
        eulerRef = 0;
601
        RmatRef = 0;
602
        Rvec = RvecOriginal;
603
        Tvec = TvecOriginal;
604
605
        determineExposure();
606
607
608
        commObj.addLog("Started reference coordinate determination.");
609
610
        CameraLibrary_EnableDevelopment();
611
        //! Initialize Camera SDK ==-
        CameraLibrary::CameraManager::X();
612
613
614
        //! At this point the Camera SDK is actively looking for all connected cameras and will initialize
615
        //! them on it's own.
616
617
        //! Get a connected camera ====
        CameraManager::X().WaitForInitialization();
618
619
        Camera *camera = CameraManager::X().GetCamera();
620
621
        //! If no device connected, pop a message box and exit ==--
622
        if (camera == 0)
623
624
            commObj.addLog("No camera found!");
625
626
627
        //! Determine camera resolution to size application window ==----
62.8
629
        int cameraWidth = camera->Width();
630
        int cameraHeight = camera->Height();
631
        camera->GetDistortionModel(distModel);
```

```
632
        cv::Mat matFrame(cv::Size(cameraWidth, cameraHeight), CV_8UC1);
633
634
        //! Set camera mode to precision mode, it directly provides marker coordinates
635
        camera->SetVideoType(Core::PrecisionMode);
636
637
        //! Start camera output ==--
638
        camera->Start();
639
640
        //! Turn on some overlay text so it's clear things are
641
        //! working even if there is nothing in the camera's view. ===---
        \ensuremath{//!} Set some other parameters as well of the camera
642
        camera->SetTextOverlay(true);
643
644
        camera->SetFrameRate(intFrameRate);
645
        camera->SetIntensity(intIntensity);
646
        camera->SetIRFilter(true);
647
        camera->SetContinuousIR(false);
648
        camera->SetHighPowerMode(false);
649
650
        //! sample some frames and calculate the position and attitude. then average those values and use that
       as zero position
651
        int numberSamples = 0;
652
        int numberToSample = 200;
        double projectionError = 0; //!< difference between the marker points as seen by the camera and the
653
       projected marker points with Rvec and Tvec
654
655
        while (numberSamples < numberToSample)</pre>
656
657
             //! Fetch a new frame from the camera ===---
658
            Frame *frame = camera->GetFrame();
659
660
            if (frame)
661
662
                 //! Ok, we've received a new frame, lets do something
663
                //! with it.
664
                if (frame->ObjectCount() == numberMarkers)
665
                     //!for(int i=0; i<frame->ObjectCount(); i++)
666
                     for (int i = 0; i < numberMarkers; i++)</pre>
667
668
                     {
669
                         cObject *obj = frame->Object(i);
670
                         list_points2dUnsorted[i] = cv::Point2d(obj->X(), obj->Y());
671
                     }
672
673
                     if (gotOrder == false)
674
                     {
675
                         determineOrder();
676
677
                    //! sort the 2d points with the correct indices as found in the preceeding order
678
       determination algorithm
679
                     for (int w = 0; w < numberMarkers; w++)</pre>
680
681
                         list_points2d[w] = list_points2dUnsorted[
      pointOrderIndices[w]];
682
683
                     list points2dOld = list points2dUnsorted;
684
685
                     //!Compute the pose from the 3D-2D corresponses
                     solvePnP(list_points3d, list_points2d,
686
      cameraMatrix, distCoeffs, Rvec, Tvec, useGuess,
      methodPNP):
687
688
                     //! project the marker 3d points with the solution into the camera image CoSy and calculate
       difference to true camera image
689
                    projectPoints(list_points3d, Rvec, Tvec,
      cameraMatrix, distCoeffs, list_points2dProjected);
690
                    projectionError = norm(list_points2dProjected,
      list points2d);
691
692
                    double maxValue = 0;
693
                     double minValue = 0;
694
                    minMaxLoc(Tvec.at<double>(2), &minValue, &maxValue);
695
                     if (maxValue > 10000 || minValue < 0)
696
697
                     {
                         ss.str("");
698
699
                         ss << "Negative z distance, thats not possible. Start the set zero routine again or
       restart Programm.";
700
                         commObj.addLog(QString::fromStdString(ss.str()));
701
                         frame->Release():
702
                         return 1;
703
                     }
704
705
                     if (projectionError > 3)
706
                         \verb|commObj.addLog("Reprojection error is bigger than 3 pixel. Correct marker| \\
707
       configuration loaded?\nMarker position measured precisely?");
```

```
frame->Release();
709
                         return 1;
710
                     }
711
                     if (norm(positionOld) - norm(Tvec) < 0.05) //!<Iterative Method needs time</pre>
712
       to converge to solution
713
714
                         add(posRef, Tvec, posRef);
715
                         add(eulerRef, Rvec, eulerRef); //!< That are not the values of yaw,</pre>
       roll and pitch yet! Rodriguez has to be called first.

numberSamples++; //!< one sample more :D
716
                         commObj.progressUpdate(numberSamples * 100 / numberToSample);
717
718
719
                     positionOld = Tvec;
720
721
                     Mat cFrame(480, 640, CV_8UC3, Scalar(0, 0, 0));
                     for (int i = 0; i < numberMarkers; i++)</pre>
722
723
724
                         circle(cFrame, Point(list_points2d[i].x,
      list_points2d[i].y), 6, Scalar(0, 225, 0), 3);
725
                     projectCoordinateFrame(cFrame);
726
72.7
                     projectPoints(list_points3d, Rvec, Tvec,
      728
729
730
                         circle(cFrame, Point(list_points2d[i].x,
      list_points2d[i].y), 3, Scalar(225, 0, 0), 3);
731
                     drawPositionText(cFrame, position,
732
      eulerAngles, projectionError);
733
734
                     QPixmap QPFrame;
735
                     QPFrame = Mat2QPixmap(cFrame);
736
                     commObj.changeImage(QPFrame);
737
                     QCoreApplication::processEvents();
738
739
740
                 frame->Release();
741
            }
742
        //! Release camera ==--
743
744
        camera->Release():
745
746
        //!Divide by the number of samples to get the mean of the reference position
747
        divide(posRef, numberToSample, posRef);
748
        divide(eulerRef, numberToSample, eulerRef); //!< eulerRef is here in Axis Angle
       notation
749
        Rodrigues(eulerRef, RmatRef);
//!-- Euler Angles, finally
750
                                                        //!< axis angle to rotation matrix
751
752
        getEulerAngles(RmatRef, eulerRef); //!< rotation matrix to euler</pre>
753
        ss.str("");
        ss << "RmatRef is:\n";
ss << RmatRef << "\n";
754
755
        ss << "Reference Position is:\n";
ss << posRef << "[mm] \n";
756
757
758
        ss << "Reference Euler Angles are:\n";
759
        ss << eulerRef << "[deg] \n";
760
761
        //! compute the difference between last obtained TVec and the average Value
762
        //! When it is large the iterative method has not converged properly so it is advised to start the
       setReference() function once again
763
        double error = norm(posRef) - norm(Tvec);
764
        if (error > 5.0)
765
            ss << "Caution, distance between reference position and last position is: " << error << "\n Start
766
       the set zero routine once again.";
767
768
        commObj.addLog(QString::fromStdString(ss.str()));
769
        commObj.progressUpdate(0);
770
        return 0;
771 }
```

Here is the call graph for this function:



3.2.2.13 startTracking()

```
int startTracking ( )
```

Start the loop that fetches frames, computes the position etc and sends it to other computers. This function is the core of this program, hence the pose estimation is done here.

Definition at line 261 of file main.cpp.

```
261
2.62
263
          gotOrder = false; //! The order of points, hence which entry in list_points3d corresponds to
264
         which in list_points2d is not calculated yet

Rvec = RvecOriginal; //! Use the value of Rvec that was set in main() as starting value
265
         for the solvePnP algorithm
266
          Tvec = TvecOriginal; //! Use the value of Tvec that was set in main() as starting value
         for the solvePnP algorithm
267
         GetLocalTime(&logDate); //! Get the current date and time to name the log file
268
269
          //! Concat the log file name as followed. The file is saved in the folder /logs in the Rigid Track
       installation folder
logFileName = "./logs/positionLog_" + QString::number(logDate.wDay) + "_" +
QString::number(logDate.wMonth) + "_" + QString::number(logDate.wYear);
logFileName += "_" + QString::number(logDate.wHour) + "_" + QString::number(logDate.wMinute) + "_" + QString::number(logDate.wSecond) + ".txt";
270
271
272
          logName = logFileName.toStdString(); //! Convert the QString to a standard string
274
         determineExposure(); //! Get the exposure where the right amount of markers is
         detected
275
276
          //! For OptiTrack Ethernet cameras, it's important to enable development mode if you
          //! want to stop execution for an extended time while debugging without disconnecting
```

```
278
        //! the Ethernet devices. Lets do that now:
279
280
        CameraLibrary_EnableDevelopment();
        CameraLibrary::CameraManager::X(); //! Initialize Camera SDK
281
282
        //! At this point the Camera SDK is actively looking for all connected cameras and will initialize
283
284
        //! them on it's own
285
286
         //! Get a connected camera
287
        CameraManager::X().WaitForInitialization();
288
        Camera *camera = CameraManager::X().GetCamera();
289
290
        //! If no camera can be found, inform user in message log and exit function
291
        if (camera == 0)
292
        {
293
             commObj.addLog("No camera found!");
294
            return 1:
295
296
297
        //! Determine camera resolution to size application window
298
        int cameraWidth = camera->Width();
299
        int cameraHeight = camera->Height();
300
        camera->SetVideoType(Core::PrecisionMode); //! Set the camera mode to precision mode, it used
301
       {\tt greyscale} \ {\tt imformation} \ {\tt for} \ {\tt marker} \ {\tt property} \ {\tt calculations}
302
303
        camera->Start(); //! Start camera output
304
305
        //! Turn on some overlay text so it's clear things are
        //! working even if there is nothing in the camera's view
306
307
        camera->SetTextOverlay(true);
308
        camera->SetExposure(intExposure);
                                                //! Set the camera exposure
        camera->SetIntensity(intIntensity); //! Set the camera infrared LED intensity camera->SetFrameRate(intFrameRate); //! Set the camera framerate to 100 Hz
309
310
311
        camera->SetIRFilter(true); //! Enable the filter that blocks visible light and only passes infrared
       liaht
        camera->SetHighPowerMode(true); //! Enable high power mode of the LEDs camera->SetContinuousIR(false); //! Disable continuous LED light
312
313
        camera->SetThreshold(intThreshold); //! Set threshold for marker detection
314
315
316
        //! Create a new matrix that stores the grayscale picture from the camera
317
        Mat matFrame = Mat::zeros(cv::Size(cameraWidth, cameraHeight), CV_8UC1);
318
        QPixmap QPFrame; //! QPixmap is the corresponding Qt class that saves images
319
        //! Matrix that stores the colored picture, hence marker points, coordinate frame and reprojected
       points
320
        Mat cFrame(480, 640, CV_8UC3, Scalar(0, 0, 0));
321
        int v = 0; //! Helper variable used to kick safety switch
322
        //! Variables for the min and max values that are needed for sanity checks
323
324
        double maxValue = 0;
325
        double minValue = 0;
326
        int framesDropped = 0; //! Ff a marker is not visible or accuracy is bad increase this counter
327
        double projectionError = 0; //! Equals the quality of the tracking
328
        setUpUDP(); //! Open sockets and ports for UDP communication
329
330
331
        if (safetyEnable) //! If the safety feature is enabled send the starting message
332
333
             //! Send enable message, hence send a 9 and then a 1
334
            data.setNum((int)(9));
335
            udpSocketSafety->write(data):
336
            data.setNum((int)(1));
337
            udpSocketSafety->write(data);
338
339
340
        //! Fetch a new frame from the camera
341
        \verb|bool| gotTime = false; //! Get the timestamp of the first frame. This time is subtracted from every |
       subseeding frame so the time starts at 0 in the logs
342
        while (!gotTime) //! While no new frame is received loop
343
        {
344
             Frame *frame = camera->GetFrame(); //! Get a new camera frame
345
             if (frame) //! There is actually a new frame
346
             {
                 timeFirstFrame = frame->TimeStamp(); //! Get the time stamp for the first frame.
347
       It is subtracted for the following frames
348
                                     //! Release the frame so the camera can continue
                 frame->Release();
349
                 gotTime = true; //! Exit the while loop
350
            }
351
        }
352
        //! Now enter the main loop that processes each frame and computes the pose, sends it and logs stuff
353
354
        while (!exitRequested) //! Check if the user has not pressed "Stop Tracking" yet
355
356
357
            Frame \starframe = camera->GetFrame(); //! Fetch a new frame from the camera
358
            if (frame) //! Did we got a new frame or does the camera still need more time
359
```


```
360
361
                framesDropped++; //! Increase by one, if everything is okay it is decreased at the end of the
       loop again
362
363
                //! Only use this frame it the right number of markers is found in the picture
364
                if (frame->ObjectCount() == numberMarkers)
365
                {
366
                     //! Get the marker points in 2D in the camera image frame and store them in the
       list_points2dUnsorted vector
367
                     //! The order of points that come from the camera corresponds to the Y coordinate
368
                     for (int i = 0; i < numberMarkers; i++)</pre>
369
370
                        cObject *obj = frame->Object(i);
371
                         list_points2dUnsorted[i] = cv::Point2d(obj->X(), obj->Y());
372
373
                    if (gotOrder == false) //! Was the order already determined? This is false for the
374
       first frame and from then on true
375
376
                        determineOrder(); //! Now compute the order
377
378
                    //! Sort the 2d points with the correct indices as found in the preceeding order
379
       determination algorithm
380
                    for (int w = 0; w < numberMarkers; w++)</pre>
381
382
                        list_points2d[w] = list_points2dUnsorted[
      pointOrderIndices[w]]; //! pointOrderIndices was calculated in determineOrder()
383
384
                     list_points2dOld = list_points2dUnsorted;
385
386
                     //! The first time the 2D-3D corresspondence was determined with gotOrder was okay.
387
                     //! But this order can change as the object moves and the marker objects appear in a
388
                     //! different order in the frame->Object() array.
389
                     //! The solution is that: When a marker point (in the camera image, hence in 2D) was at
390
                     //! a position then it wont move that much from one frame to the other.
391
                     //! So for the new frame we take a marker object and check which marker was closest this
       point
392
                     //! in the old image frame? This is probably the same (true) marker. And we do that for
       every other marker as well.
393
                     //! When tracking is good and no frames are dropped because of missing markers this should
       work every frame.
394
                    for (int j = 0; j < numberMarkers; j++)</pre>
395
                     {
396
                        minPointDistance = 5000; //! The sum of point distances is set to
       something unrealistic large
397
                        for (int k = 0; k < numberMarkers; k++)</pre>
398
399
                             //! Calculate N_2 norm of unsorted points minus old points
                             currentPointDistance = norm(
400
      list_points2dUnsorted[pointOrderIndices[j]] -
      list_points2dOld[k]);
401
                             //! If the norm is smaller than minPointDistance the correspondence is more likely
       to be correct
                             if (currentPointDistance <</pre>
402
      minPointDistance)
403
404
                                 //! Update the array that saves the new point order
405
                                minPointDistance =
      currentPointDistance:
406
                                 pointOrderIndicesNew[j] = k;
407
                             }
408
                        }
409
                    }
410
411
                     //! Now the new order is found, set the point order to the new value
412
                     for (int k = 0; k < numberMarkers; k++)</pre>
413
414
                        pointOrderIndices[k] = pointOrderIndicesNew[k];
415
                         list_points2d[k] = list_points2dUnsorted[
      pointOrderIndices[k]];
416
417
                     //! Save the unsorted position of the marker points for the next loop
418
                    list points2dOld = list points2dUnsorted;
419
420
421
                     //!Compute the object pose from the 3D-2D corresponses
422
                     solvePnP(list_points3d, list_points2d,
      cameraMatrix, distCoeffs, Rvec, Tvec, useGuess,
      methodPNP):
423
424
                    //! Project the marker 3d points with the solution into the camera image CoSy and calculate
       difference to true camera image
425
                    projectPoints(list_points3d, Rvec, Tvec,
      cameraMatrix, distCoeffs, list_points2dProjected);
426
                    projectionError = norm(list_points2dProjected,
      list_points2d); //! Difference of true pose and found pose
```

```
428
                                           //! Increase the framesDropped variable if accuracy of tracking is too bad
429
                                                (projectionError > 5)
430
431
                                                   framesDropped++;
432
433
                                          else
434
                                          {
435
                                                   framesDropped = 0; //! Set number of subsequent frames dropped to zero because error
               is small enough and no marker was missing
436
437
438
                                           //! Get the min and max values from TVec for sanity check
439
                                          minMaxLoc(Tvec.at < double > (2), & minValue, & maxValue);
440
441
                                           //! Sanity check of values. negative z means the marker CoSy is behind the camera, that's
              not possible.
442
                                           if (minValue < 0)</pre>
443
444
                                                   commObj.addLog("Negative z distance, that is not possible. Start the set zero
               routine again or restart Program.");
445
                                                   frame->Release(); //! Release the frame so the camera can move on
                                                   camera->Release(); //! Release the camera
446
                                                  {\tt closeUDP} (); //! Close all UDP connections so the programm can be closed later
447
              on and no resources are locked
448
                                                 return 1; //! Exit the function
449
450
451
                                           //! Next step is the transformation from camera CoSy to navigation CoSy
                                          //! Compute the relative object position from the reference position to the current one //! given in the camera CoSy: f\ T_C^{NM} = Tvec - Tvec_{Ref} \f
452
453
454
                                          subtract(Tvec, posRef, position);
455
456
                                           //! Transform the position from the camera CoSy to the navigation CoSy with INS alligned
              457
458
            position:
459
                                          position = V;
                                                                             //! Position is the result of the preceeding calculation
                                          position[2] *= invertZ; //! Invert Z if check box in GUI is activated,
460
              hence height above ground is considered
461
                                           //! Realtive angle between reference orientation and current orientation
462
463
                                          Rodrigues (Rvec, Rmat); //! Convert axis angle respresentation to ordinary rotation
              matrix
464
465
                                           //! The difference of the reference rotation and the current rotation
466
                                           //! \f$ R_{ NM } = M_{ NC } \times R_{ CM } \f$
467
                                          Rmat = RmatRef.t() *Rmat;
468
469
                                           //! Euler Angles, finally
                                           getEulerAngles(Rmat, eulerAngles); //! Get the euler angles
470
               from the rotation matrix % \left( \frac{1}{2}\right) =\frac{1}{2}\left( \frac{1}{2}\right) =\frac{1}{2}\left
471
                                          eulerAngles[2] += headingOffset; //! Add the heading offset to the
              heading angle
472
                                           //! Compute the velocity with finite differences. Only use is the log file. It is done here
               because the more precise time stamp can be used
474
                                          frameTime = frame->TimeStamp() - timeOld; //! Time between the old frame
               and the current frame
475
                                          timeOld = frame->TimeStamp();
                                                                                                               //! Set the old frame time to the current one
                                          velocity[0] = (position[0] - positionOld[0]) /
476
             frameTime; //! Calculate the x velocity with finite differences
                                           velocity[1] = (position[1] - positionOld[1]) /
477
             frameTime; //! Calculate the y velocity with finite differences
478
                                          velocity[2] = (position[2] - positionOld[2]) /
             479
              next frame velocity calcuation
480
481
                                           //! Send position and Euler angles over WiFi with 100 Hz
482
                                           sendDataUDP(position, eulerAngles);
483
                                           //! Save the values in a log file, values are:
484
                                           //! Time sinc tracking started Position
485
                                                                                                                                      Euler Angles
                                                                                                                                                                        Velocity
                                           logfile.open(logName, std::ios::app); //! Open the log file, the folder is
486
               RigidTrackInstallationFolder/logs
            logfile << frame->TimeStamp() - timeFirstFrame << ";" << position[0] << ";" << position[2] << ";";" << position[2] << ";";" << position[2] << ";";" << eulerAngles[0] << ";" << eulerAngles[1] << ";" << eulerAngles[2] << ";";";";
487
488
                                          logfile << velocity[0] << ";" << velocity[1] << ";" <<
489
             velocity[2] << "\n";</pre>
490
                                          logfile.close(); //! Close the file to save values
491
492
                                 //! Check if the position and euler angles are below the allowed value, if wes send OKAY signal
493
```

```
(1), if not send shutdown signal (0)
494
                //! Absolute x, y and z position in navigation CoSy must be smaller than the allowed distance
495
                if (safetyEnable)
496
497
                    if ((abs(position[0]) < safetyBoxLength && abs(position[1]) <</pre>
      safetyBoxLength && abs(position[2]) < safetyBoxLength))</pre>
498
499
                        //! Absolute Euler angles must be smaller than allowed value. Heading is not considered
500
                        if ((abs(eulerAngles[0]) < safetyAngle && abs(eulerAngles[1]) <</pre>
      safetyAngle))
501
502
                            //! Send the OKAY signal to the desired computer every 5th time
503
504
                                data.setNum((int)(1));
505
                                udpSocketSafety\hbox{->write(data); //! Send the 1}
506
                                v = 0; //! reset the counter that is needed for decimation to every 5th time
       step
507
508
                        //! The euler angles of the object exceeded the allowed euler angles, send the shutdown
509
       signal (0)
510
                        else
511
                            data.setNum((int)(0)); //! Send the shutdown signal, a 0
512
513
                            udpSocketSafety->write(data);
514
                            commObj.addLog("Object exceeded allowed Euler angles, shutdown signal sent."
      ); //! Inform the user
515
516
517
518
                    //! The position of the object exceeded the allowed position, shut the object down
519
520
521
                        data.setNum((int)(0)); //! Send the shutdown signal, a 0
                        udpSocketSafety->write(data);
522
                        commObj.addLog("Object left allowed area, shutdown signal sent."); //! Inform
523
       the user
524
525
526
                }
527
                //! Inform the user if tracking system is disturbed (marker lost or so) or error was too big
528
529
                if (framesDropped > 10)
530
531
                    if (safetyEnable) //! Also send the shutdown signal
532
                        data.setNum((int)(0)); //! Send the shutdown signal, a 0
533
                        udpSocketSafety->write(data);
534
535
536
                    commObj.addLog("Lost marker points or precision was bad!"); //! Inform the user
537
                    framesDropped = 0;
538
539
                //! Rasterize the frame so it can be shown in the GUI
540
                frame->Rasterize(cameraWidth, cameraHeight, matFrame.step,
541
      BACKBUFFER_BITSPERPIXEL, matFrame.data);
542
543
                //! Convert the frame from greyscale as it comes from the camera to rgb color
544
                cvtColor(matFrame, cFrame, COLOR_GRAY2RGB);
545
546
                //! Project (draw) the marker CoSy origin into 2D and save it in the cFrame image
547
                projectCoordinateFrame(cFrame);
548
549
                //! Project the marker points from 3D to the camera image frame (2d) with the computed pose
550
                projectPoints(list_points3d, Rvec, Tvec,
      551
552
553
                    //! Draw a circle around the projected points so the result can be better compared to the
       real marker position
554
                    //! In the resulting picture those are the red dots
555
                    circle(cFrame, Point(list_points2d[i].x,
      list_points2d[i].y), 3, Scalar(225, 0, 0), 3);
556
               }
557
558
                //! Write the current position, attitude and error values as text in the frame
559
                drawPositionText(cFrame, position, eulerAngles, projectionError);
560
                //! Send the new camera picture to the GUI and call the GUI processing routine
561
                QPixmap QPFrame;
562
563
                QPFrame = Mat2QPixmap(cFrame);
                commObj.changeImage(QPFrame); //! Update the picture in the GUI
564
565
                QCoreApplication::processEvents(); //! Give Qt time to handle everything
566
                //! Release the camera frame to fetch the new one
567
                frame->Release();
568
```

```
569 }
570 }
571
572 //! User choose to stop the tracking, clean things up
573 closeUDP(); //! Close the UDP connections so resources are deallocated
574 camera->Release(); //! Release camera
575 return 0;
576 |
```

Here is the call graph for this function:



Here is the caller graph for this function:



3.2.2.14 testAlgorithms()

```
void testAlgorithms ( )
```


Project some points from 3D to 2D and then check the accuracy of the algorithms. Mainly to generate something that can be shown in the camera view so the user knows everything loaded correctly.

Definition at line 952 of file main.cpp.

```
953 {
954
       int _methodPNP;
956
957
       std::vector<Point2d> noise(numberMarkers);
958
959
       RvecOriginal = Rvec;
960
       TvecOriginal = Tvec;
961
962
       projectPoints(list_points3d, Rvec, Tvec, cameraMatrix,
     distCoeffs, list_points2dProjected);
963
964
       ss.str("");
965
       ss << "Unsorted Points 2D Projected \n";
966
       ss << list_points2dProjected << "\n";
967
       commObj.addLog(QString::fromStdString(ss.str()));
968
       Mat cFrame(480, 640, CV_8UC3, Scalar(0, 0, 0));
for (int i = 0; i < numberMarkers; i++)</pre>
969
970
971
972
            circle(cFrame, Point(list_points2dProjected[i].x, list_points2dProjected[i].y), 6, Scalar(0, 255, 0
     ), 3);
973
974
975
       projectCoordinateFrame(cFrame);
976
       ss.str("");
978
        979
980
       ss << list_points2dProjected << "\n";
981
982
       randn(noise, 0, 0.5);
       add(list_points2dProjected, noise, list_points2dProjected);
983
984
985
       ss << "======== With Noise Points =======\n";
986
        ss << list_points2dProjected << "\n";
987
       commObj.addLog(QString::fromStdString(ss.str()));
988
989
990
       bool useGuess = true;
       _methodPNP = 0; //! < 0 = iterative 1 = EPNP 2 = P3P 4 = UPNP //! < not used
991
992
993
        solvePnP(list_points3d, list_points2dProjected, cameraMatrix,
     distCoeffs, Rvec, Tvec, useGuess, _methodPNP);
994
995
       ss.str("");
996
       ss << "====
997
       ss << "=======\n";
       ss << "rvec: " << "\n";
998
       ss << Rvec << "\n";
ss << "tvec: " << "\n";
999
1000
        ss << Tvec << "\n";
1001
1002
1003
        commObj.addLog(QString::fromStdString(ss.str()));
1004
         _methodPNP = 1; //!< 0 = iterative 1 = EPNP 2 = P3P 4 = UPNP UPnP not used
1005
        Rvec = cv::Mat::zeros(3, 1, CV_64F);
1006
     Tvec = cv::Mat::zeros(3, 1, cv_64F);
solvePnP(list_points3d, list_points2dProjected, cameraMatrix,
distCoeffs, Rvec, Tvec, useGuess, _methodPNP);
1007
1008
1009
        ss.str("");
1010
        ss << "-----\n":
1011
        ss << "======\n";
1012
        ss << "rvec: " << "\n";
1013
        ss << Rvec << "\n";
ss << "tvec: " << "\n";
1014
1015
        ss << Tvec << "\n";
1016
1017
1018
        projectPoints(list_points3d, Rvec, Tvec, cameraMatrix,
     distCoeffs, list_points2dProjected);
       for (int i = 0; i < numberMarkers; i++)</pre>
```

```
circle(cFrame, Point(list_points2dProjected[i].x, list_points2dProjected[i].y), 3, Scalar(255, 0, 0
      ), 3);
1022
        QPixmap QPFrame;
1023
1024
         QPFrame = Mat2QPixmap(cFrame);
         commObj.changeImage(QPFrame);
1025
1026
         QCoreApplication::processEvents();
1027
        commObj.addLog(QString::fromStdString(ss.str()));
1028
        if (numberMarkers == 4)
1029
              methodPNP = 2; //! < 0 = iterative 1 = EPNP 2 = P3P 4 = UPNP //! < not used
1030
            Rvec = cv::Mat::zeros(3, 1, CV_64F);
Tvec = cv::Mat::zeros(3, 1, CV_64F);
1031
1032
1033
             solvePnP(list_points3d, list_points2dProjected,
      cameraMatrix, distCoeffs, Rvec, Tvec, useGuess, _methodPNP);
1034
1035
             ss.str("");
            ss << "======\n";
1036
            ss << "======\n";
1037
1038
             ss << "rvec: " << "\n";
            ss << Rvec: " << "\n";
ss << "tvec: " << "\n";
1039
1040
            ss << Tvec << "\n";
1041
1042
            projectPoints(list_points3d, Rvec, Tvec, cameraMatrix,
1043
     distCoeffs, list_points2dProjected);
1044
           for (int i = 0; i < numberMarkers; i++)</pre>
1045
             {
1046
                circle(cFrame, Point(list_points2dProjected[i].x, list_points2dProjected[i].y), 3, Scalar(255,
      0, 0), 3);
1047
1048
             double projectionError = norm(list_points2dProjected, list_points2d);
1049
             putText(cFrame, "Testing Algorithms Finished", cv::Point(5, 420), 1, 1, cv::Scalar(255, 255, 255));
             drawPositionText(cFrame, position, eulerAngles, projectionError)
1050
1051
1052
            QPixmap QPFrame;
1053
            QPFrame = Mat2QPixmap(cFrame);
1054
             commObj.changeImage(QPFrame);
1055
             QCoreApplication::processEvents();
1056
            commObj.addLog(QString::fromStdString(ss.str()));
1057
        }
1058
1059
         _methodPNP = 4; //! < 0 = iterative 1 = EPNP 2 = P3P 4 = UPNP //! < not used
1060
         Rvec = cv::Mat::zeros(3, 1, CV_64F);
1061
         Tvec = cv::Mat::zeros(3, 1, CV_64F);
         solvePnP(list_points3d, list_points2dProjected, cameraMatrix,
1062
      distCoeffs, Rvec, Tvec, useGuess, _methodPNP);
1063
1064
         ss.str("");
1065
         ss << "======\n";
         ss << "========
                                             =======\n";
1066
                                      UPNP
         ss << "rvec: " << "\n";
1067
        ss << Rvec: " << "\n";
ss << Rvec << "\n";
ss << "tvec: " << "\n";
1068
1069
1070
        ss << Tvec << "\n";
1071
1072
         commObj.addLog(QString::fromStdString(ss.str()));
1073
1074
         Rvec = RvecOriginal:
         Tvec = TvecOriginal;
1075
1076
1077 }
```

Here is the call graph for this function:



Here is the caller graph for this function:



3.2.3 Variable Documentation

3.2.3.1 commObj

commObject commObj

class that handles the communication from main.cpp to the GUI

Now declare variables that are used across the main.cpp file. Basically almost every variable used is declared here.

Definition at line 68 of file main.cpp.

3.3 RigidTrack/RigidTrack.cpp File Reference

Rigid Track GUI source that contains functions for GUI events.

```
#include "RigidTrack.h"
#include <QProcess>
#include <QdesktopServices>
#include <QDir>
#include <QWessageBox>
#include 'QUrl>
#include "main.h"
#include "communication.h"
#include dependency graph for RigidTrack.cpp:
```



3.3.1 Detailed Description

Rigid Track GUI source that contains functions for GUI events.

Author

Florian J.T. Wachter

Version

1.0

Date

April, 8th 2017


3.4 RigidTrack/RigidTrack.h File Reference

Rigid Track GUI source header with Qt Signals and Slots.

```
#include <QtWidgets/QMainWindow>
#include "ui_RigidTrack.h"
#include <qpixmap.h>
#include "main.h"
#include "communication.h"
Include dependency graph for RigidTrack.h:
```



This graph shows which files directly or indirectly include this file:



3.4.1 Detailed Description

Rigid Track GUI source header with Qt Signals and Slots.

Author

Florian J.T. Wachter

Version

1.0

Date

April, 8th 2017

Index

| BACKBUFFER_BITSPERPIXEL | main.h, 60 |
|-----------------------------|-------------------------|
| main.cpp, 40 | l A l |
| 1 B 10 B 11 | eulerAngles |
| calcBoardCornerPositions | main.cpp, 42 |
| main.cpp, 15 | eulerRef |
| calibrateCamera | main.cpp, 42 |
| main.cpp, 16 | exitRequested |
| main.h, 54 | main.cpp, 43 |
| calibrateGround | fromoTimo |
| main.cpp, 17 | frameTime |
| main.h, 55 | main.cpp, 43 |
| camera_started | GET OPTIMIZED |
| main.cpp, 40 | precomp.hpp, 77 |
| cameraMatrix | getEulerAngles |
| main.cpp, 41 | • |
| closeUDP | main.cpp, 23 |
| main.cpp, 19 | gotOrder |
| main.h, 57 | main.cpp, 43 |
| commObj | headingOffset |
| main.cpp, 41 | _ |
| main.h, 74 | main.cpp, 43 |
| coordinateFrame | IDI ICON1 |
| main.cpp, 41 | resource.h, 78 |
| coordinateFrameProjected | IPAdressObject |
| main.cpp, 41 | main.cpp, 44 |
| currentMinIndex | main.h, 75 |
| main.cpp, 41 | IPAdressSafety |
| currentPointDistance | main.cpp, 44 |
| main.cpp, 41 | main.h, 75 |
| тат.орр, т | |
| data | IPAdressSafety2 |
| main.cpp, 42 | main.cpp, 44 |
| datagram | main.h, 75 |
| main.cpp, 42 | intExposure |
| Debug/moc_RigidTrack.cpp | main.cpp, 43 |
| QT_MOC_LITERAL, 7 | intFrameRate |
| Debug/moc_communication.cpp | main.cpp, 43 |
| QT MOC LITERAL, 5 | intIntensity |
| | main.cpp, 44 |
| determineExposure | intThreshold |
| main.cpp, 20 | main.cpp, 44 |
| main.h, 57 | invertZ |
| determineOrder | main.cpp, 44 |
| main.cpp, 22 | main.h, 75 |
| main.h, 59 | |
| distCoeffs | list_points2d |
| main.cpp, 42 | main.cpp, 45 |
| distModel | list_points2dDifference |
| main.cpp, 42 | main.cpp, 45 |
| drawPositionText | list_points2dOld |
| main.cpp, 23 | main.cpp, 45 |

74 INDEX

| list_points2dProjected | IPAdressSafety2, 44 |
|------------------------------|-----------------------------|
| main.cpp, 45 | intExposure, 43 |
| list_points2dUnsorted | intFrameRate, 43 |
| main.cpp, 45 | intIntensity, 44 |
| list_points3d | intThreshold, 44 |
| main.cpp, 45 | invertZ, 44 |
| loadCalibration | list_points2d, 45 |
| main.cpp, 24 | list points2dDifference, 45 |
| main.h, 61 | list points2dOld, 45 |
| loadCameraPosition | list_points2dProjected, 45 |
| main.cpp, 25 | list points2dUnsorted, 45 |
| main.h, 61 | list_points3d, 45 |
| loadMarkerConfig | loadCalibration, 24 |
| main.cpp, 25 | loadCameraPosition, 25 |
| main.h, 62 | loadMarkerConfig, 25 |
| logDate | logDate, 46 |
| main.cpp, 46 | |
| logFileName | logFileName, 46 |
| main.cpp, 46 | logName, 46 |
| logName | logfile, 46 |
| | M_CN, 46 |
| main.cpp, 46 logfile | M_HeadingOffset, 46 |
| main.cpp, 46 | main, 27 |
| тат. с рр, 40 | Mat2QPixmap, 28 |
| M CN | methodPNP, 47 |
| main.cpp, 46 | minPointDistance, 47 |
| M_HeadingOffset | numberMarkers, 47 |
| main.cpp, 46 | pointOrderIndices, 47 |
| main | pointOrderIndicesNew, 47 |
| main.cpp, 27 | portObject, 47 |
| main.cpp | portSafety, 48 |
| BACKBUFFER BITSPERPIXEL, 40 | portSafety2, 48 |
| - | posRef, 48 |
| calcBoardCornerPositions, 15 | position, 48 |
| calibrateCamera, 16 | positionOld, 48 |
| calibrateGround, 17 | projectCoordinateFrame, 29 |
| camera_started, 40 | Rmat, 48 |
| cameraMatrix, 41 | RmatRef, 49 |
| closeUDP, 19 | Rvec, 49 |
| commObj, 41 | RvecOriginal, 49 |
| coordinateFrame, 41 | safety2Enable, 49 |
| coordinateFrameProjected, 41 | safetyAngle, 49 |
| currentMinIndex, 41 | safetyBoxLength, 49 |
| currentPointDistance, 41 | , , |
| data, 42 | safetyEnable, 50 |
| datagram, 42 | sendDataUDP, 29 |
| determineExposure, 20 | setHeadingOffset, 30 |
| determineOrder, 22 | setReference, 30 |
| distCoeffs, 42 | setUpUDP, 32 |
| distModel, 42 | ss, 50 |
| drawPositionText, 23 | startStopCamera, 33 |
| eulerAngles, 42 | startTracking, 34 |
| eulerRef, 42 | strBuf, 50 |
| exitRequested, 43 | testAlgorithms, 39 |
| frameTime, 43 | timeFirstFrame, 50 |
| getEulerAngles, 23 | timeOld, 50 |
| gotOrder, 43 | Tvec, 50 |
| headingOffset, 43 | TvecOriginal, 51 |
| IPAdressObject, 44 | udpSocketObject, 51 |
| IPAdressSafety, 44 | udpSocketSafety, 51 |
| | |

INDEX 75

| udpSocketSafety2, 51 | main.cpp, 48 |
|--------------------------------|--|
| useGuess, 51 | position |
| velocity, 51 | main.cpp, 48 |
| main.h | positionOld |
| calibrateCamera, 54 | main.cpp, 48 |
| calibrateGround, 55 | precomp.hpp |
| closeUDP, 57 | GET_OPTIMIZED, 77 |
| commObj, 74 | projectCoordinateFrame |
| determineExposure, 57 | main.cpp, 29 |
| determineOrder, 59 | main.h, 63 |
| drawPositionText, 60 | |
| IPAdressObject, 75 | qCleanupResources_RigidTrack |
| IPAdressSafety, 75 | qrc_RigidTrack.cpp, 9 |
| IPAdressSafety2, 75 | qInitResources_RigidTrack |
| invertZ, 75 | qrc_RigidTrack.cpp, 9 |
| loadCalibration, 61 | QT_MOC_LITERAL |
| loadCameraPosition, 61 | Debug/moc_RigidTrack.cpp, 7 |
| loadMarkerConfig, 62 | Debug/moc_communication.cpp, 5 |
| methodPNP, 75 | Release/moc_RigidTrack.cpp, 8 |
| portObject, 75 | Release/moc_communication.cpp, 6 |
| portSafety, 76 | QT_RCC_MANGLE_NAMESPACE |
| portSafety2, 76 | qrc_RigidTrack.cpp, 8 |
| projectCoordinateFrame, 63 | QT_RCC_PREPEND_NAMESPACE |
| safety2Enable, 76 | qrc_RigidTrack.cpp, 8 |
| safetyAngle, 76 | qrc_RigidTrack.cpp |
| safetyBoxLength, 76 | qCleanupResources_RigidTrack, 9 |
| safetyEnable, 76 | qInitResources_RigidTrack, 9 |
| sendDataUDP, 64 | QT_RCC_MANGLE_NAMESPACE, 8 |
| setHeadingOffset, 64 | QT_RCC_PREPEND_NAMESPACE, 8 |
| setReference, 65 | Release/moc_RigidTrack.cpp |
| setUpUDP, 67 | QT_MOC_LITERAL, 8 |
| startStopCamera, 67 | Release/moc_communication.cpp |
| startTracking, 68 | QT_MOC_LITERAL, 6 |
| testAlgorithms, 74 | resource.h |
| Mat2QPixmap | IDI ICON1, 78 |
| main.cpp, 28 | RigidTrack/ modelest.h, 3 |
| methodPNP | RigidTrack/GeneratedFiles/Debug/moc_RigidTrack.cpp, |
| main.cpp, 47 | 7 |
| main.h, 75 minPointDistance | RigidTrack/GeneratedFiles/Debug/moc_communication. |
| | cpp, 5 |
| main.cpp, 47 | RigidTrack/GeneratedFiles/Release/moc RigidTrack.← |
| numberMarkers | cpp, 7 |
| main.cpp, 47 | RigidTrack/GeneratedFiles/Release/moc_communication. |
| | cpp, 6 |
| pointOrderIndices | RigidTrack/GeneratedFiles/qrc_RigidTrack.cpp, 8 |
| main.cpp, 47 | RigidTrack/GeneratedFiles/ui_RigidTrack.h, 10 |
| pointOrderIndicesNew | RigidTrack/RigidTrack.cpp, 78 |
| main.cpp, 47 | RigidTrack/RigidTrack.h, 79 |
| portObject | RigidTrack/communication.cpp, 3 |
| main.cpp, 47 | RigidTrack/communication.h, 4 |
| main.h, 75 | RigidTrack/main.cpp, 11 |
| portSafety | RigidTrack/main.h, 52 |
| main.cpp, 48 | RigidTrack/precomp.hpp, 77 |
| main.h, 76 | RigidTrack/resource.h, 78 |
| portSafety2 | Rmat |
| main.cpp, 48 | main.cpp, 48 |
| main.h, 76 | RmatRef |
| posRef | main.cpp, 49 |

76 INDEX

Rvec main.cpp, 51 main.cpp, 49 velocity RvecOriginal main.cpp, 51 main.cpp, 49 safety2Enable main.cpp, 49 main.h, 76 safetyAngle main.cpp, 49 main.h, 76 safetyBoxLength main.cpp, 49 main.h, 76 safetyEnable main.cpp, 50 main.h, 76 sendDataUDP main.cpp, 29 main.h, 64 setHeadingOffset main.cpp, 30 main.h, 64 setReference main.cpp, 30 main.h, 65 setUpUDP main.cpp, 32 main.h, 67 main.cpp, 50 startStopCamera main.cpp, 33 main.h, 67 startTracking main.cpp, 34 main.h, 68 strBuf main.cpp, 50 testAlgorithms main.cpp, 39 main.h, 74 timeFirstFrame main.cpp, 50 timeOld main.cpp, 50 Tvec main.cpp, 50 TvecOriginal main.cpp, 51 udpSocketObject main.cpp, 51 udpSocketSafety main.cpp, 51 udpSocketSafety2 main.cpp, 51

useGuess